

SPIKED DRINK TESTING

Rapid Response™ 3-in-1 Spiked Drink Drug Test Kit (GHB, Ketamine and Benzo)

| For beverage safety testing use only.



In social settings where beverage tampering or spiking is a concern, testing your beverage is a defensive step to protect against potential harm.

The Rapid Response™ 3-in-1 Spiked Drink Drug Test Kit is an advanced multi-drug test capable of detecting three common drugs used for drug-facilitated sexual assault (date-rape drugs) in beverages: Benzodiazepines at 2µg/mL, Ketamine 10µg/mL, and GHB 500µg/mL.

Details

Sample: Beverage
Format: Cassette
Time-to-Result: 5 minutes
Storage Condition: 35.6-86°F / 2-30°C
Test Principle: Lateral Flow Assay



Detects 3 Common Drugs

Benzodiazepines, ketamine, and GHB are commonly associated with drug-facilitated sexual assault and can be unknowingly slipped into drinks.



Easy & Convenient

A portable test includes a pipette for sample collection and a buffer tube, for a mess-free and easy testing process.



Kit Content

Each kit includes 10 test cassettes individually packaged with droppers, 10 sample dilution buffer tubes, a quick reference guide, and a product insert.

What are Benzodiazepines?

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called γ aminobutyric acid (GABA). At higher doses, benzodiazepines can produce dissociation effects and anterograde amnesia during the effects. Benzodiazepines also have little to no taste or odor.

What is Ketamine?







Ketamine is a dissociative anesthetic developed in 1963 to replace PCP (Phencyclidine). While Ketamine is still used in human anesthesia and veterinary medicine, it is becoming increasingly abused as a street drug. Ketamine is molecularly similar to PCP and thus creates similar effects including numbness, loss of coordination, sense of invulnerability, muscle rigidity, aggressive / violent behavior, slurred or blocked speech, exaggerated sense of strength, and a blank stare. Ketamine is colorless and odorless.

What is GHB?

γ -Hydroxybutyric acid (GHB) is an endogenous metabolite in the brain and peripheral organs. It has many characteristics of a neurotransmitter and has been studied for potential therapeutic use in the treatment of narcolepsy, drug addiction, and symptoms of withdrawal and to induce anesthesia. At higher doses, GHB produces sedation and a trance-like state with loss of memory. GHB has little smell or taste.

3-in-1 Spiked Drink Drug Test Kit


Test Procedure

<p>1 Remove the test cassette from the sealed pouch.</p> 	<p>2 Peel open the aluminum foil on the buffer tube.</p> 	<p>3 Use the dropper to draw a tube of beverage. Transfer 2 drops of beverage (approx. 80µL) to the buffer tube.</p> 
<p>4 Close the lid of the buffer tube.</p> 	<p>5 Shake the buffer tube from side to side 5-8 times.</p> 	<p>6 Place the test cassette on a clean and level surface. Transfer 3 drops (approx. 120µL) to the specimen well (S) of the test cassette. Start the timer. Read the results at 5 minutes.</p> 


These instructions are for illustrative purposes only. Please read the instructions supplied with the test before use.

Result Interpretation

Negative



Positive



Invalid

For BZO/KET: No line appears in the Control region (C).

For GHB: If the color pad has a purple color before applying the beverage sample, do not use the test.

Ordering Information

Product Code	Product Name	Contents
D3.39-8C29-10	Rapid Response™ 3-in-1 Spiked Drink Drug Test Kit	10 Tests/Case

To learn more contact your local **Sales Representative**, call us at **+1 888 506-2658**, or email us at **info@lochnessmedical.com**