REF VITD-13C1

## Rapid Response®

### Vitamin D Test Cassette

(Fingerstick Whole Blood)

#### Product Insert

For self-testing *in vitro* diagnostic use.

A rapid test for the semi-quantitative detection of 25-hydroxyvitamin D in human fingerstick whole blood.

#### Intended Use

The Rapid Response® Vitamin D Test Cassette is a rapid chromatographic immunoassay for the semi-quantitative detection of 25-hydroxyvitamin D (25 (OH) D) in human fingerstick whole blood. This assay provides a preliminary diagnostic test result and can be used to screening for Vitamin D deficiency.

#### Summary

Vitamin D refers to a group of fat-soluble secosteroids responsible for increasing intestinal absorption of calcium, iron, magnesium, phosphate and zinc. In humans, the most important compounds in this group are vitamin D3 and vitamin D2. Vitamin D3 is naturally produced in the human skin through the exposure to ultraviolet light and Vitamin D2 is mainly obtained from foods. Vitamin D is transported to the liver where it is metabolized to 25-hydroxy Vitamin D. In medicine, a 25-hydroxy Vitamin D blood test is used to determine Vitamin D concentration in the body. The blood concentration of 25-hydroxy Vitamin D (including D2 and D3) is considered the best indicator of Vitamin D status. Vitamin D deficiency is now recognized as a global epidemic.<sup>2</sup> Virtually every cell in our body has Receptors for Vitamin D, meaning that they all require "Sufficient" Level of Vitamin D for adequate functioning. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin deficiency has been linked to various serious diseases: Osteoporosis, Osteomalacia, Multiple Sclerosis, Cardiovascular Diseases, Pregnancy Complications, Diabetes, Depression, Strokes, Autoimmune Diseases, Flu, Different Cancers, Infectious Diseases, Alzheimer, Obesity and Higher Mortality etc.<sup>3</sup>

#### Principle

The Rapid Response® Vitamin D Test Cassette is an immunoassay based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip. During testing, 25 (OH) D present in the specimen will compete with 25 (OH) D on the test line for limited amount of anti-25 OH Vitamin D antibodies in the conjugate. The higher the concentration of 25 (OH) D in the specimen, the lighter the T line would be. The result should be read according to the Colour Card provided with the kit. To serve as a procedural control, a colour line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

#### Precautions

Please read all the information in this product insert before performing the test.

- For self-testing *in vitro* diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with doctor or medical professional.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle or touch the test window of the test cassette.
- The kit must not be frozen or used after the expiration date printed on the package.
- Keep out of the reach of children.
- The used test should be discarded according to local regulations.

#### Storage and Stability

Store as packaged in the sealed pouch either at room temperature or refrigerated 2-30°C (36-86°F). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

#### Limitations

- The Rapid Response® Vitamin D Test Cassette provides only a semi-quantitative analytical result. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- Other clinically available tests are required if questionable results are obtained.

#### Quality Control

A procedural control is included in the test. A coloured line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct procedural technique.

#### Performance Characteristics

##### Accuracy

A clinical evaluation was conducted comparing the results obtained using the Rapid Response® Vitamin D Test Cassette to predicate Device (Vitamin D Rapid Test). The in-house clinical trial included 90 whole blood specimens. The results demonstrated with an overall accuracy of 94.4%.

Method	Predicate Device (Vitamin D Test)			Total Results
	Deficient	Insufficient	Sufficient	
Rapid Response® Vitamin D Test Cassette	Deficient	4	3	7
	Insufficient	0	53	55
	Sufficient	0	0	28
<b>Total Results</b>	<b>4</b>	<b>56</b>	<b>30</b>	<b>90</b>
<b>Accuracy</b>	<b>&gt;99.9%</b>	<b>94.6%</b>	<b>93.3%</b>	<b>94.4%</b>

#### Extra Informations

##### 1. How does the Vitamin D test work?

In medicine, a 25-hydroxy Vitamin D is the main storage form of vitamin D in the body. Therefore, the overall status of vitamin D can be determined by detecting the content of 25-hydroxy Vitamin D. 25-hydroxy Vitamin D level less than 30 ng/mL in case of a positive result, indicates Vitamin D Deficiency or Insufficiency. Vitamin D supplements can be recommended in these cases.

##### 2. When should the test be used?

The clinical application of 25-hydroxy Vitamin D is mainly for diagnosis, treatment and monitoring of rickets (children), osteomalacia, postmenopausal osteoporosis and renal osteopathy. Vitamin D deficiency is also associated with many other diseases, including cancer, cardiovascular disease, autoimmune diseases, diabetes and depression. Monitor your vitamin D levels to determine whether to take vitamin D supplements. The Rapid Response® Vitamin D Test can be used any time of the day.

##### 3. Can the result be incorrect?

The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the Rapid Response® Vitamin D Test Cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Due to the immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always recommended for such tests based on immunological principles.

## Rapid Response®

### Vitamin D Test Cassette

(Fingerstick Whole Blood)

#### Setting Up the Test

Wash your hands with soap and rinse with clear warm water.

Unpack the test components from the kit and make sure that all the packaging is intact.


For each test you will need:


Test Cassette

Product Insert  
(This Document)


Buffer

Lancet

Alcohol Pad

Capillary Dropper

Timer  
(Not provided)


Colour Card

#### Procedure


Bring the pouch to room temperature before opening it. Open the pouch, remove the test cassette and place it on a clean and level surface. Run the test within one hour and best results will be obtained if the test is performed immediately after opening the foil pouch.


Remove the dropper, buffer vial, lancet and Alcohol pad, place them close to the test cassette.


Carefully pull off and dispose the released cap of the lancet.


Use the provided Alcohol pad to clean the fingertip of the middle or ring finger as the puncture site. Allow to air dry.


Press the lancet, on the side from where the cap was extracted; the tip retracts automatically and safely after use. Without touching the puncture site, massage the hand towards the fingertip of the middle or ring finger to be punctured.

Keeping the hand down massage the end that was pricked to obtain a blood drop.


Without squeezing the capillary dropper bulb, put it in contact with the blood. **The blood migrates into the capillary dropper through the capillarity to the line indicated on the capillary dropper.**


You may massage again your finger to obtain more blood if the blood does not reach the indicated line. Avoid of air bubbles.


**Release the blood collected into the Specimen well (S) of the cassette, by squeezing the dropper bulb.**

Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and **add 2 drops of buffer into the Buffer well (B) of the cassette and start a timer.**


Wait for the coloured line(s) to appear. **Read results at 10 minutes.** Compare the T line intensity with "Vitamin D Colour Card" provided with the kit to get the Vitamin D level in your blood. Do not interpret the result after 20 minutes.

#### Reading the Results

Please refer to the illustration and compare the T line intensity with "Vitamin D Colour Card" provided with the kit.

**Note: Do not mix use Vitamin D colour card from different batches.**

##### 4. How to interpret the test if the colour and the intensity of the lines are different?

Please refer to the illustration and compare the T line intensity with "Vitamin D Colour Card" provided with the kit.

##### 5. If I read the result after 20 minutes, will the result be reliable?

No. The result should be read at 10 minutes after adding the buffer. The result is unreliable after 20 minutes.

##### 6. What do I have to do if the result is deficient or insufficient?

If the result is deficient or insufficient, it means that the Vitamin D level in blood is less than 30 ng/mL and that you should consult a physician. The physician will decide whether additional analysis should be performed.


##### 7. What do I have to do if the result is sufficient?

If the result is sufficient, it means that the Vitamin D level is higher than or equal to 30 ng/mL and is within the normal range. A case of Vitamin D toxicity (hypercalcemia), though rare, but cannot be excluded based on such test results. However, if the symptoms persist, it is recommended to consult a physician.

#### Bibliography

- Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health. *Mayo Clinic Proceedings*. 81 (3): 353-73.
- Eriksen EF, Gluerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. *Biogerontology*. 3 (1-2): 73-7.
- Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. *Alternative Medicine Review*. 10 (2): 94-111.

#### Glossary of Symbols

**BTNX, Inc.**  
722 Rosebank Road  
Pickering, ON, L1W 4B2, Canada  
Technical Support: 1-888-339-9964

25-OH Vitamin D Level	Deficient
Reference Range (ng/mL)	0-10
Reference Range (nmol/L)	0-25


**DEFICIENT: Two coloured lines appear.** One is in the control region (C) and another should be in the test region (T). The intensity of line in the test region (T) is **equal to Deficient line (0-10 ng/mL) on the colour card**, it indicates Vitamin D level is deficient.

25-OH Vitamin D Level	Insufficient
Reference Range (ng/mL)	10-30
Reference Range (nmol/L)	25-75


**INSUFFICIENT: Two coloured lines appear.** One is in the control region (C) and another should be in the test region (T). The intensity of line in the test region (T) is **equal to Insufficient line (10-30 ng/mL) on the colour card**, it indicates Vitamin D level is insufficient.

25-OH Vitamin D Level	Sufficient
Reference Range (ng/mL)	30-100
Reference Range (nmol/L)	75-250

**SUFFICIENT: Two coloured lines appear.** One line should be always in the control region (C) and faint coloured line appears in the test region (T). The intensity of line in the test region (T) is **equal to Sufficient line (30-100 ng/mL) on the colour card**, it indicates Vitamin D level is sufficient.

**EXCESS: One coloured line appears in the control line region (C). No coloured line appears in the test line region (T).** It indicates Vitamin D level maybe excess, it is recommended to consult a physician.

**INVALID: Control line fails to appear.** Incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor

**BTNX Inc.**  
Tel: +1-888-339-9964  
Fax: +1-905-944-0406  
Email/courriel: support@btnx.com  
[www.btnx.com](http://www.btnx.com)

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#### Références

1. Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health. *Mayo Clinic Proceedings*. 81 (3): 353-73.
2. Eriksen EF, Gluerp H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. *Biogerontology*. 3 (1-2): 73-7.
3. Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. *Alternative Medicine Review*. 10 (2): 94-111.

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**BTNX, Inc.**  
722 Rosebank Road  
Pickering, ON, L1W 4B2, Canada  
Support Technique: 1-888-339-9964

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