Rapid Response®

Rapid Response® 4-in-1 COVID-19, FLU A/B & RSV Test -

At Home

REF CFR-19PHTC1 REF CFR-19PHTC4
REF CFR-19PHTC2 REF CFR-19PHTC5

Quick Reference Instructions


A rapid test for the detection of SARS-CoV-2, Influenza A/B and RSV antigens in anterior nasal specimens.

For *in vitro* diagnostic use only.

Carefully read the instructions before performing the test. Failure to follow the instructions may result in inaccurate test results.

Contents

NOTE:

Do not use this test on children under 2 years of age. Children between 2 and 18 years of age must be aided or supervised by an adult when carrying out the test.


Prepare for the Test

Bring the device(s), reagents and specimens to room temperature (15~30°C; 59~86°F) before use.

- 1 Wash your hands with soap and water or use hand sanitizer for 20 seconds.
- 2 Check the expiration date on the box. Do not use the kit if it has been damaged or has expired.
- 3 Remove the test device from its packaging. For the best results, the assay should be performed within one hour.

Test Procedure

- 1 Remove the test pen and the base from the packaging.
USE THE TEST WITHIN 1 HOUR OF OPENING.
- 2 Remove the swab cap from the test pen.
Discard the Swab Cap.
- 3 Gently insert the swab end of the pen in one nostril about 1/2 - 3/4 of an inch.
Firmly rub the swab at least **5 times** against the inside walls of the nostril in a circular motion. Do not just spin the swab.
Make sure you have removed the swab cap (see step 2).


- 4 Repeat the process with the same swab in the other nostril.

Did you swab both nostrils? If not, inaccurate results can occur.

Test Procedure

- 5a Place the base on a flat surface. Place the swab end of the pen into the base.
Locate white band on the pen.
- 5b Press firmly on the test pen to insert it all the way into the base so that the white band is completely covered by the base.
White band covered.

WARNING: Failure to insert the test pen all the way into the base can lead to inaccurate results.

- 6 Leaving the test upright, set a timer and read the results at 15 minutes.
WARNING: Read the results at 15 minutes. Do not read the results after 30 minutes.
After the test is completed, dispose of used materials in the trash.

Follow the advice of your local provincial or territory health department for guidance or confirmation testing if necessary, and if unwell seek medical assistance.

Read and Interpret Your Results

For COVID-19 Test

- POSITIVE (+)**
- Two coloured lines appear on the membrane. One line appears in the control region (C) and another line appears in the test region (T).
- NOTE:** A positive test result means it is very likely you currently have COVID-19 disease. Please contact your provincial or territory health department for guidance on virus testing services to get a laboratory PCR test immediately.

Read and Interpret Your Results

For COVID-19 Test

- NEGATIVE (-)**
- Only one coloured line appears, in the control region (C). No apparent coloured line appears in the test region (T).
- NOTE:** A negative result for COVID-19 does not mean you do not have COVID-19. If you have symptoms, you should follow the guidance from your local provincial or territory health departments, and if unwell seek medical assistance. Negative results are presumptive and may need to be confirmed with a molecular assay. If symptoms continue or infection is suspected, you should be tested again with at least 24 hours and no more than 48 hours between tests as SARS-CoV-2 antigens cannot be precisely detected in all phases of an infection.

- INVALID**
- No coloured line appears in the control region (C), whether a test line(s) is present or not. Results from any test which has not produced a control line at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

For Influenza A/B Test

- INFLUENZA A POSITIVE (+)**
- One coloured line appears in the control region (C), and another coloured line appears in the A region (A).
- INFLUENZA B POSITIVE (+)**
- One coloured line appears in the control region (C), and another coloured line appears in the B region (B).

For Influenza A/B Test


INFLUENZA A+B POSITIVE (+)

One coloured line appears in the control region (C), and two other coloured lines appear in both the A region (A) and the B region (B).

Co-infection with influenza A and B is rare. A clinical specimen that generates positive results for both A and B should be considered an invalid result, and another test should be performed. If the test is again positive for both influenza A and B, the specimen should be re-tested by another method prior to reporting of results.

A positive test result means that the virus that causes influenza A or influenza B was detected in your sample, and it is very likely that you have influenza. Please continue to observe local hygiene and safety measures.

NOTE: There is a very small chance that this test can give a result that is incorrect (a false positive). A positive result does not rule out co-infections with other pathogens or identify any specific influenza virus subtype.

NEGATIVE (-)

Only one coloured line appears in the control region (C), and a line appears neither in the A region (A) nor B region (B).

A negative test result means it is unlikely you have influenza A/B disease. Please continue to observe local hygiene and safety measures.

NOTE: A negative result does not mean you do not have influenza, and if symptoms persist, you should seek medical attention and further testing if required.

INVALID

No coloured line appears in the control region (C), whether a test line(s) is present or not.

Results from any test which has not produced a control line at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

For RSV Test

POSITIVE (+)

Two coloured lines appear on the membrane. One line appears in the control region (C) and another line appears in the test region (T).

A positive test result means that the virus that causes RSV was detected in your sample, and it is very likely that you have RSV. Please continue to observe local hygiene and safety measures.

NOTE: There is a very small chance that this test can give a result that is incorrect (a false positive). A positive result does not rule out co-infections with other pathogens or identify any specific RSV subtype.

NEGATIVE (-)

Only one coloured line appears, in the control region (C). No apparent coloured line appears in the test region (T).

A negative test result means it is unlikely you have RSV disease. Please continue to observe local hygiene and safety measures.

NOTE: A negative result does not mean you are not infectious or do not have RSV. If symptoms persist, you should seek medical attention and further testing if required.

INVALID

No coloured line appears in the control region (C), whether a test line(s) is present or not. Results from any test which has not produced a control line at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

Intended Use

The Rapid Response® 4-in-1 COVID-19, FLU A/B & RSV Test – At Home for self-testing is an *in vitro* immunoassay. The assay is intended for home testing (or self-testing). Children aged between 2 and 18 years old must be supervised or tested by an adult when carrying out the test. The assay is an *in vitro* immunochromatographic assay for the qualitative detection of antigens of SARS-CoV-2, Influenza A virus, Influenza B virus and Respiratory Syncytial Virus (RSV) in nasal swab specimens collected from individuals with the respiratory infection for COVID-19, influenza A/B, and/or RSV within the first 7 days of the onset of symptoms. The assay obtains a preliminary result only aiding in the diagnosis of COVID-19, Influenza A (FLU A), Influenza B (FLU B), and/or RSV. This test has not been cleared for use in asymptomatic individuals. Testing results should not be the sole basis for treatment or other management decisions.

Limitations of the Test

1. The Rapid Response® 4-in-1 COVID-19, FLU A/B & RSV Test – At Home is for self-testing *in vitro* diagnostic use, and should only be used for the qualitative detection of viral antigens specific for SARS-CoV-2, Influenza A virus, Influenza B virus, and RSV. The intensity of colour in a positive line should not be evaluated as "quantitative or semi-quantitative".

2. Both viable and nonviable viruses are detectable with the kit.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
4. Failure to follow the TEST PROCEDURE and RESULT INTERPRETATION may adversely affect test performance and/or invalidate the test result.
5. Negative results may not mean you are not infectious and if symptoms are present, you must seek immediate further testing.
6. A negative result does not rule out infection with another type of respiratory virus.
7. Repeat testing (e.g. within 1-3 days) is recommended if infection is still suspected, in a high-risk setting, or for occupational or other requirements.
8. The test is less reliable in the later phase of infection (more than 7 days after the onset of COVID-19, Influenza A/B and RSV symptoms) and in asymptomatic individuals.
9. A positive result cannot necessarily determine whether you are infectious.
10. Caution should be taken when inserting the swab into the nasal cavity.

Storage and Stability

- Store the test kit at 2-30°C (35.6-86°F) when not in use.
- **DO NOT FREEZE.**
- Kit contents are stable until the expiration dates marked on the outer packaging and container.
- Do not use after the expiry date on the packaging.
- Do not use the kit when a component is missing or damaged.

Safety Information

- Wearing a mask can help protect you and those around you if you are in an area with community transmission and physical distancing is not possible.
- Follow the directions of your local provincial or territory government health department to protect yourself.
- Test kit solutions should only be used as directed; do not ingest; do not dip the swab into provided solution or other liquid before inserting the swab into the nose; avoid contact with skin and eyes; keep out of the reach of children and pets before and after use.
- If the extraction buffer comes in contact with the skin or eyes, flush with plenty of water. If irritation persists, seek medical advice from a doctor or your local medical centre.

Frequently Asked Questions

Q. What are COVID-19, Influenza A/B and RSV?

- A.** COVID-19 is an acute respiratory infectious disease caused by the SARS-CoV-2 virus, a novel beta coronavirus. SARS-CoV-2 mostly spreads person-to-person, both by individuals with symptoms of COVID-19 and by infected people without symptoms. Based on current knowledge, the incubation period is 1-4 days, mostly 4-5 days. The symptoms include fever, fatigue, and cough. Influenza A is a highly contagious respiratory disease caused by RNA viruses belonging to the Myxoviridae family. It is transmitted by droplets and comes in contact with contaminated objects. The incubation period is 1-4 days and the infection period can be up to 5-7 days after symptom onset. Typical symptoms include sudden high fever (>38°C; >100.4°F), muscle pain, headache, dry cough, and sore throat. In severe cases, complications, such as pneumonia and myocarditis, may develop and even lead to death. Influenza B belongs to the Bunyaviridae family and usually presents with an acute onset, with chills and fever, and a rapid rise in body temperature within a few hours to 24 hours. RSV belongs to the pneumovirus family and is a single negative-stranded RNA virus transmitted by close contact or droplets, with an incubation period of 2-8 days. In infants, the symptoms are bronchiolitis (wheezing, breathing difficulties, and feeding difficulties), and in

adults, the symptoms of the common cold are common. Severe infections can cause pneumonia, and are the leading cause of hospitalization in infants.

Q. What does an invalid test result mean?

- A.** If no control line appears on the test, the result is invalid (even if any test line appears). An invalid result means that the test was not able to tell if you had COVID-19, Influenza A/B, and RSV. If the test is invalid, a new swab should be used to collect a new nasal specimen and the test should be repeated using all new test components.

Q. Will this test hurt?

- A.** No, the nasal swab is not sharp and should not cause injuries. The swab may feel slightly uncomfortable. If you feel pain, please stop the test and seek advice from your healthcare provider.

Q. What are the known potential risks and benefits of this test?

- A.** Potential risks include:
- Possible discomfort during sample collection.
 - Possible incorrect test results (see the Results Interpretation section).
- Potential benefits include:
- The results, along with other information, can help you and your healthcare provider make informed decisions about your care.
 - The results of this test may help limit the spread of COVID-19, Influenza A/B and RSV to your family and others in your community.

Q. What is the difference between an antigen test and a molecular test?

- A.** There are different types of tests for COVID-19, Influenza A/B, and RSV. Molecular tests (also known as PCR tests) detect genetic material from the virus. Antigen tests, such as the Rapid Response® 4-in-1 COVID-19, FLU A/B & RSV Test – At Home detect proteins from the virus. Antigen tests are very specific for SARS-CoV-2, Influenza A virus, Influenza B virus and RSV, but are not as sensitive as molecular tests. This means that a positive result is highly accurate, whereas a negative result does not rule out infection.

Q. How accurate is this test?

- A.** The performance of the Rapid Response® 4-in-1 COVID-19, FLU A/B & RSV Test – At Home was established in a prospective clinical study using an EUA authorized molecular test as a comparator method (PPA (90.0%) and NPA (99.6%) for COVID-19, PPA (92.1%) and NPA (99.7%) for Influenza A, PPA (92.0%) and NPA (99.9%) for Influenza B, PPA (90.9%) and NPA (99.9%) for RSV).

Q. What if you test positive?


- A.** A positive test result for COVID-19 means that antigens from the virus that cause COVID-19 were detected, and it is very likely you currently have COVID-19. A positive test result for Influenza A means that antigens from the virus that cause Influenza A were detected, and it is very likely that you currently have Influenza A. A positive test result for Influenza B indicates that antigens from the virus that causes Influenza B were detected, and it is very likely you currently have Influenza B. A positive test result for RSV indicates that antigens from the virus that causes RSV were detected, and it is very likely that you currently have RSV. There is a very small chance that this test can give a positive result that is wrong (a false positive result). If you test positive, you should self-isolate at home per CDC recommendations to stop the virus from spreading to others. Seek follow-up care with healthcare providers immediately. Your healthcare provider will work with you to determine how best to care for you based on your test result(s) along with your medical history, and your symptoms.

Q. What if you test negative?

- A.** All COVID-19 and/or Influenza A/B and/or RSV antigen test negative results are presumptive, and confirmation with a molecular assay may be necessary. If you continue to have

symptoms of COVID-19 and/or Influenza A/B and/or RSV and both the first and second tests are negative, you may not have COVID-19, Influenza A/B, and/or RSV. However, you should follow up with your healthcare provider. Symptomatic individuals who test negative should repeat testing at least twice over three days with at least 48 hours between tests and at least three times over five days with at least 48 hours between tests if they are asymptomatic.

Glossary of Symbols

Scan

for more detailed instructions found in the complete Instructions for Use or visit:
www.btnx.com/product/covid-fluab-rsv-test-at-home

BTNX, Inc.
722 Rosebank Road
Pickering, ON L1W 4B2 Canada
Technical Support: 1-888-339-9964
Email: support@btnx.com
Mon.-Fri. 9:00AM-5:30PM ET
www.btnx.com

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