Rapid Response™

COVID-19 Antigen Rapid Test Device

REF:COV-19C25

Product Insert (EN)

Intended Use

The Rapid Response™ COVID-19 Antigen Rapid Test is an *in vitro* immunochromatographic assay for the direct and qualitative detection of SARS-CoV-2 viral nucleoprotein antigens from human nasopharyngeal secretions, nasal and oropharyngeal secretions from individuals suspected of COVID-19. The test is intended for use with Nasopharyngeal or Oropharyngeal swab specimens within the first two weeks of symptom onset, and with Nasal swab specimens within 7 days from symptom onset. This test is intended for professional use only.

Principle

The Rapid Response™ COVID-19 Antigen Rapid Test detects SARS-CoV-2 viral antigens through visual interpretation of color development. Anti-SARS-CoV-2 antibodies are immobilized on the test region of the nitrocellulose membrane. Anti-SARS-CoV-2 antibodies conjugated to colored particles are immobilized on the conjugated pad. A sample is added to the extraction buffer which is optimized to release the SARS-CoV-2 antigens from specimen. During testing, the extracted antigens bind to anti-SARS-CoV-2 antibodies conjugated to colored particles. As the specimen migrates along the strip by capillary action and interacts with reagents on the membrane, the complex will be captured by the anti-SARS-CoV-2 antibodies at the test region. Excess colored particles are captured at the internal control zone. The presence of a colored band in the test region indicates a positive result for the SARS-CoV-2 viral antigens, while its absence indicates a negative result. A colored band at the control region serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking is working.

Precautions

- For *in vitro* Diagnostic Use Only.
- Read the Package Insert prior to use. Directions should be read and followed carefully.
- Do not use kit or components beyond the expiration date.
- The device contains material of animal origin and should be handled as a potential biohazard. Do not use if pouch is damaged or open.
- Test devices are packaged in foil pouches that exclude moisture during storage. Inspect each foil pouch before opening. Do not use devices that have holes in the foil or where the pouch has not been completely sealed. Erroneous result may occur if test reagents or components are improperly stored. Do not use the Extraction Buffer if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection during specimen collection and assaying.
- All patient specimens should be handled and discarded as if they are biologically hazardous. All specimens must be mixed thoroughly before testing to ensure a representative sample prior to testing.
- All specimens must be mixed thoroughly before testing to ensure a representative sample prior to testing.
- Failure to bring specimens and reagents to room temperature before testing may decrease assay sensitivity. Inaccurate or inappropriate specimen collection, storage, and transport may yield false negative test results.
- Do not add specimen to the reaction area (result area).
- To avoid contamination, do not touch the reaction area (result area).

- Use a new extraction tube and nozzle for each specimen to avoid cross-contamination.
- Avoid skin contact with buffer.
- If infection with SARS-CoV-2 is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions and sent to state or local health departments for testing.
- Viral isolation in cell culture and initial characterization of viral agents recovered in cultures of SARS-CoV-2 specimens are NOT recommended, except in a BSL3 laboratory using BSL3 work practices.

Materials

Materials provided

- Individually packed test devices
- Extraction buffer
- Extraction tube
- Nozzle with filter
- Individually packed swabs
- Tube stand
- Product insert

Materials required but not provided

- Clock, timer, or stopwatch

Storage and Stability

- The Rapid Response™ COVID-19 Antigen Rapid Test Device kit should be stored at 35.6~86°F (2~30°C) when not in use.
- **DO NOT FREEZE.**
- Test contents are stable until the expiration dates marked on the foil pouch and containers.
- Do not use kit or components beyond the expiration date.
- Test cassettes are packaged in foil pouches that exclude moisture during storage. Inspect each foil pouch before opening. Do not test kit components if the primary packaging is damaged or if the pouch has not been completely sealed. Erroneous result may occur if test reagents or components are improperly stored. Protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination can lead to inaccurate results.

Collection of Specimens

Nasopharyngeal swab (NP swab):

1. Remove the swab from its packaging.
2. Gently insert the sterile swab into the nostril parallel to the palate, not upward. The distance should be equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx, or until resistance is encountered. Gently rub and roll the swab, rotating against the nasal wall. Leave in place several seconds to saturate tip with secretions. Slowly remove the swab while rotating it.
3. Process the swab as soon as possible after collecting the specimen.

Oropharyngeal swab (OP swab):

1. Remove the swab from its packaging.
2. Insert the swab completely from the mouth into the throat, centering on the red part of the throat wall and maxillary tonsils and rub the bilateral throat tonsils and throat wall moderately. Avoid touching the tongue and remove the swab.
3. Process the swab as soon as possible after collecting the specimen.

Nasal swab (N swab):

1. Remove the swab from its packaging.
2. Tilt patient's head back 70°. Insert the swab through the anterior nares in contact with nasal septum at least 0.5 inches inside the nostril until mild resistance is encountered at the middle turbinate.


3. Using a circular motion, the nasal orifice should be swabbed for a minimum of five seconds.
4. Compress the nostril with the fingers to trap the swab tip and rotate the tip for a minimum of five seconds.
5. Remove and repeat for the other nostril with the same swab.
6. Process the swab as soon as possible after collecting the specimen.

Note:

1. Use only synthetic fiber swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit further testing.
2. Swabs specimens should be tested as soon as possible after collection. Use freshly collected specimens for best test performance.
3. If not tested immediately, swab specimens may be stored at 35.6-46.4°F (2-8°C) for 24 hours after collection.
4. Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of the sample with the interpretation of test results.

Test Procedure

Bring devices, reagents, and specimens and/or controls to room temperature (59~86°F; 15~30°C) before use.

1. Label a clean extraction tube with patient or control identification. And place it into the tube stand.
2. Gently mix extraction buffer. Without touching the buffer bottle to the extraction tube, add 10 drops into the extraction tube. *OR* Empty the contents of one extraction buffer vial into the collection tube.
3. Insert the swab with the collected specimen into the extraction tube. Swirl the swab, mixing well. Squeeze the swab 10-15 times by compressing the walls of the tube against the swab.
4. **Let the solution stand for 2 minutes.**
5. Remove the swab while pressing the swab head firmly against the inner wall of the tube to release as much liquid as possible. Dispose of the used swab in accordance with the appropriate biohazard waste disposal protocol.
6. For each specimen, open the foil pouch just before testing and remove the test device and put it on a clean, level surface. For best results, the assay should be performed within one hour. Label the test device with patient or control identification.
7. Attach nozzle to sample extraction tube. Invert the tube and add 3 drops of the extracted solution into the sample well of the test device by gently squeezing the tube.
8. Start the timer. Wait for coloured line(s) to appear. Read results at 15 minutes.

Results Interpretation

POSITIVE: Two coloured lines appear on the membrane. One line appears in the control region (C) and another line appears in the test region (T)*.

NEGATIVE: Only one coloured line appears, in the control region (C). No apparent coloured line appears in the test region (T).

INVALID: Control line fails to appear. Results from any test which has not produced a control line at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

NOTE:

1. The colour intensity in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of colour in the test region should be considered positive. Note that this is a qualitative test only and cannot determine the concentration of analytes in the specimen.
2. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

Quality Control

Internal Procedural Controls

The Rapid Response™ COVID-19 Antigen Rapid Test has built-in (procedural) controls. Each test device has an internal standard zone to ensure proper sample flow. The user should confirm that the coloured line located at the "C" region is present before reading the result.

External Positive and Negative Controls

Good laboratory practice suggests testing positive and negative external controls to ensure that the test reagents are working and that the test is correctly performed.

Limitations

1. The Rapid Response™ COVID-19 Antigen Rapid Test is for professional *in vitro* diagnostic use only. It should only be used for the qualitative detection of SARS-CoV-2 viral nucleoprotein antigens in human nasal, nasopharyngeal and oropharyngeal specimen.
2. This is a qualitative test. The intensity of colour in a positive line should not be evaluated as "quantitative or semi-quantitative". A quantitative value or rate of increase/decrease in the concentration of SARS-CoV-2 viral nucleoprotein antigens cannot be determined with this test.
3. Both viable and nonviable SARS-CoV-2 viruses are detectable with Rapid Response™ COVID-19 Antigen Rapid Test.
4. The Rapid Response™ COVID-19 Antigen Rapid Test Device detects the presence of SARS-CoV-2 viral nucleoprotein antigens in specimen and should not be used as the sole criterion for a diagnosis of COVID-19. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, rather, the results of this test should be interpreted in conjunction with all other clinical and laboratory findings available to the physician.
5. Operators must closely follow the instructions in the "Specimen Collection and Storage", "Test Procedure" and "Result Interpretation" sections. Failure to follow the instructions may adversely affect test performance or lead to inaccurate test results.
6. Results obtained with this assay, particularly in the case of weak test lines that are difficult to interpret, should be used in conjunction with other clinical information available to the physician.
7. Negative results do not preclude SARS-CoV-2 infection and should be

8. confirmed via molecular assay. The performance of the device has not been assessed on specimens from individuals who have been infected with emerging variants of SARS-CoV-2 of public health concern.
9. Analytical studies with recombinant protein shows that this device could detect the recombinant nucleocapsid protein of SARS-CoV-2 variants of B.1.1.7, B.1.351, P1 and B.1.617.2.

Performance Characteristics

Analytical Sensitivity (Limit of Detection):

The limit of detection was determined with a quantified SARS-CoV-2 virus and has been evaluated at $2 \times 10^{2.4}$ TCID₅₀/mL. The limit of detection was also determined with recombinant SARS-CoV-2 nucleoprotein and has been evaluated at 0.4 ng/mL.

Clinical Evaluation:

Study 1: With Nasal Swab as a sample type:

Clinical performance characteristics for Rapid Response™ COVID-19 Antigen Rapid Test Device were evaluated in a multi-site prospective, single-blind, controlled clinical trial in the USA from October 2020 to December 2020. The study was performed by Point-of-Care operators with no laboratory experience. These clinical trials were aimed to evaluate the performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device with a Nasal Swab by comparing with an FDA EUA approved RT-PCR comparator.

The performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device was established with direct anterior nasal swabs collected from individual symptomatic patients (within 7 days from onset) who were suspected of COVID-19. Two samples from each patient were collected – one for PCR and another for Rapid Antigen test. 51 positive specimens and 128 negative specimens were confirmed by RT-PCR.

Table 1: Rapid Response™ COVID-19 Antigen Rapid Test Clinical Evaluation with Nasal Swabs:

		RT-PCR		Total
		Positive	Negative	
Rapid Response™ COVID-19 Antigen Rapid Test Device	Positive	46	0	46
	Negative	5	128	133
Total		51	128	179

Diagnostic Sensitivity: 90.2% (78.6% ~ 96.2%)*
 Diagnostic Specificity: 100.0% (96.5% ~ 100.0%)*
 Overall Agreement: 97.2% (93.2% ~ 98.9%)*
 *95% Confidence Interval

Study 2: With Nasopharyngeal Swab as a sample type:

Clinical trials for Rapid Response™ COVID-19 Antigen Rapid Test Device were performed at two Point-of-care sites in the USA from August 2020 to October 2020. These clinical trials were aimed to evaluate the performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device by comparing with an FDA EUA approved RT-PCR comparator.

The performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device was established with 82 direct nasopharyngeal swabs collected and enrolled from individual symptomatic patients (within 14 days) who were suspected of COVID-19. Samples were freshly collected from 2 sites where the operators were minimally trained. 46 positive specimens and 36 negative specimens were confirmed by RT-PCR.

The performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device based on the results from these two sites is summarized below.

Table 2: Site 1 for Nasopharyngeal Swab Specimen vs. RT-PCR

		RT-PCR		Total
		Positive	Negative	

Rapid Response™ COVID-19 Antigen Rapid Test Device	Positive	23	0	23
	Negative	1	21	22
Total		24	21	45

Diagnostic Sensitivity: 95.8% (79.8% ~ 99.3%)*
 Diagnostic Specificity: 100.0% (84.5% ~ 100.0%)*
 Overall Agreement: 97.8% (88.4% ~ 99.6%)*
 *95% Confidence Interval

Table 3: Site 2 for Nasopharyngeal Swab Specimen vs. RT-PCR

		RT-PCR		Total
		Positive	Negative	
Rapid Response™ COVID-19 Antigen Rapid Test Device	Positive	21	0	21
	Negative	1	15	16
Total		22	15	37

Diagnostic Sensitivity: 95.5% (78.2% ~ 100.0%)*
 Diagnostic Specificity: 100.0% (79.6% ~ 100.0%)*
 Overall Agreement: 100.0% (86.2% ~ 99.5%)*
 *95% Confidence Interval

Table 4: Site 1 & 2 Combined for Nasopharyngeal Swab Specimen vs. RT-PCR

		RT-PCR		Total
		Positive	Negative	
Rapid Response™ COVID-19 Antigen Rapid Test Device	Positive	44	0	44
	Negative	2	36	38
Total		46	36	82

Diagnostic Sensitivity: 95.6% (85.5% ~ 98.8%)*
 Diagnostic Specificity: 100.0% (90.4% ~ 100.0%)*
 Overall Agreement: 97.6% (91.5% ~ 99.3%)*
 *95% Confidence Interval

Study 3: With Nasopharyngeal or Oropharyngeal Swab as sample type:

Clinical evaluation was performed to compare the results obtained by Rapid Response™ COVID-19 Antigen Rapid Test Device and a RT-PCR comparator assay. Patients who presented with symptoms within 14 days were included in the study.

The performance was evaluated for 85 prospective clinical specimens which included 43 specimens of nasopharyngeal swab and 42 specimens of oropharyngeal swab. 55 positive specimens and 30 negative specimens were confirmed by RT-PCR.

Another 131 retrospective clinical specimen who are either asymptomatic or with mild symptoms (within 14 days of symptom onset) were collected to confirm the specificity. All of these are confirmed negative through PCR. 59 specimens were of nasopharyngeal swab and 72 specimens of oropharyngeal swab.

The performance of Rapid Response™ COVID-19 Antigen Rapid Test Device results based on the various parameters is summarized below:

Table 5: Oropharyngeal Swab Specimen vs. RT-PCR Positive

Days from onset of symptoms	PCR positive at any time	Rapid Response™ COVID-19 Antigen Rapid Test Device		
		Antigen Positive	PPA	95% Confidence Interval
≤7	6	6	100.00%	51.68% to 100%
8-14	24	22	91.67%	71.52% to 98.54%
≤14	30	28	93.33%	76.49% to 98.83%

Table 6: Nasopharyngeal Swab Specimen vs. RT-PCR Positive

Days from onset of symptoms	PCR positive at any time	Rapid Response™ COVID-19 Antigen Rapid Test Device		
		Antigen Positive	PPA	95% Confidence Interval
≤7	4	4	100.00%	39.57% to 100%
8-14	21	20	95.24%	74.13% to 99.75%
≤14	25	24	96.00%	77.67% to 99.79%

Table 7: Oropharyngeal and Nasopharyngeal Swab Specimen vs. RT-PCR Positive

Days from onset of symptoms	PCR positive at any time	Rapid Response™ COVID-19 Antigen Rapid Test Device		
		Antigen Positive	PPA	95% Confidence Interval
≤7	10	10	100%	65.55% to 100%
8-14	45	42	93.33%	80.69% to 98.26%
≤14	55	52	94.55%	83.93% to 98.58%

Table 8: Summary of Rapid Response™ COVID-19 Antigen Rapid Test Clinical Evaluation:

		RT-PCR		Total
		Positive	Negative	
Rapid Response™ COVID-19 Antigen Rapid Test Device	Positive	52	0	52
	Negative	3	161	164
Total		55	161	216

Relative Sensitivity: 94.55% (83.93%~98.58%)
 *Relative Specificity: 100% (97.1%~100%)
 Overall Agreement: 98.61% (95.66%~99.64%)
 *95% Confidence Interval

Cross Reactivity:

Cross reactivity with the following organisms has been studied. Samples positive for the following organisms were found negative when tested with the Rapid Response™ COVID-19 Antigen Rapid Test (Nasopharyngeal/Oropharyngeal Swab).

HCoV-HKU1	Norovirus
HCoV-OC43	Influenza A (H5N1)
HCoV-NL63	Influenza A (H7N9)
HCoV-229E	Influenza A (H7N7)
Measles virus	Influenza B Victoria lineage
Streptococcus pneumoniae	Influenza B Yamaqata lineage
Epstein-Barr virus	Respiratory syncytial virus
Coxsackie virus A16	Adenovirus
Bordetella Para pertussis	Legionella pneumophila
Influenza A (H1N1) pdm09	Mycoplasma pneumoniae
Influenza A (H3N2)	Chlamydia pneumoniae
Parainfluenza 1/2/3 virus	Streptococcus pyogenes
Human metapneumovirus	Streptococcus agalactiae
Rhinovirus	Group C Streptococcus
Mump virus	Staphylococcus aureus

Interfering Substances

The following substances, naturally present in respiratory specimens or that may be artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. None of them were found to affect test performance of the Rapid Response™ COVID-19 Antigen Rapid Test Device.

3 OTC nasal sprays: 10%	Phenylephrine: 10 mg/ml
3 OTC mouthwashes: 10%	Phenylpropanolamine: 20 mg/ml
3 OTC throat drops: 10%	Relenza® (zanamivir): 20 mg/ml
4-acetamidophenol: 10 mg/ml	Rimantadine: 500 ng/ml
Acetylsalicylic acid: 20 mg/ml	Tamiflu® (oseltamivir): 100 mg/ml
Albuterol: 20 mg/ml	Tobramycin: 40 mg/ml

Dextromethorphan: 10 mg/ml	Chlorpheniramine: 5 mg/ml
Diphenhydramine: 5 mg/ml	Dexamethasone: 5 mg/ml
Doxylamine succinate: 1 mg/ml	Mupirocin: 250 µg/ml
Flunisolide: 3 mg/ml	Oxymetazoline: 10 mg/ml
Guaiacol glyceryl ether: 20 mg/ml	Triamcinolone: 14 mg/ml
Mucin: 1%	

Bibliography

- Forni, D., Cagliani, R., Clerici, M. & Sironi, M. Molecular evolution of human coronavirus genomes. Trends Microbiol. 25, 35–48 (2017).
- Ithete, N. L. et al. Close relative of human Middle East respiratory syndrome coronavirus in bat, South Africa. Emerg. Infect. Dis. 19, 1697–1699 (2013).

Glossary of Symbols

Consult instructions for use Test per Kit Authorized Representative

Store between 35.6°F to 86°F Use by Do Not Reuse

Lot Number For *in vitro* diagnostic use only Catalogue #

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BTNX Inc.
 722 Rosebank Road,
 Pickering, ON L1W 4B2
 Canada

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MDSS GmbH