

Rapid Response™

Filariasis IgG/IgM Test Cassette (Whole Blood/Serum/Plasma)

REF FIL-13C40

Product Insert

A rapid test for the qualitative detection of IgG and IgM antibodies to Filariasis parasites (W. Bancrofti and B. Malayi) in human whole blood, serum, or plasma.

For laboratory *in vitro* diagnostic use only.

Intended Use

The Rapid Response™ Filariasis IgG/IgM Test Cassette is a rapid chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to Filariasis parasites (W. Bancrofti and B. Malayi) in whole blood, serum, or plasma to aid in the diagnosis of Filariasis infection.

Summary

The lymphatic filariasis known as Elephantiasis, mainly caused by W. bancrofti and B. malayi, affects about 120 million people over 80 countries^{1,2}. The disease is transmitted to humans by the bites of infected mosquitoes within which the microfilariae sucked from an infected human subject develops into third-stage larvae. Generally, repeated and prolonged exposure to infected larvae is required for establishment of human infection. The definitive parasitologic diagnosis is the demonstration of microfilariae in blood samples³. However, this gold standard test is restricted by the requirement for nocturnal blood collection and lack of adequate sensitivity. Detection of circulating antigens is commercially available. Its usefulness is limited for W. bancrofti⁴. In addition, microfilaremia and antigenemia develop from months to years after exposure.

Antibody detection provides an early means to detect filarial parasite infection. Presence of IgM to the parasite antigens suggest current infection, whereas, IgG corresponds to late stage of infection or past infection⁵. Furthermore, identification of conserved antigens allows 'pan-filaria' test to be applicable. Utilization of recombinant proteins eliminates cross-reaction with individuals having other parasitic diseases⁶. The Rapid Response™ Filariasis IgG/IgM Test Cassette uses conserved recombinant antigens to simultaneously detect antibody to the W. bancrofti and B. malayi parasites without the restriction on specimen collection.

Principle

The Rapid Response™ Filariasis IgG/IgM Test Cassette is a qualitative membrane-based immunoassay for the detection of IgG and IgM antibodies to Filariasis in whole blood, serum, or plasma. In this test procedure, anti-human IgG is immobilized

in the IgG line region of the test, anti-human IgM is immobilized in the IgM line region of the test. After specimen is added to the specimen well of the device, it reacts with Filariasis antigen coated particles in the test. This mixture migrates chromatographically along the length of the test and interacts with the immobilized anti-human IgG/ anti-human IgM in respective line. If the specimen contains Filariasis antibodies, colored line will appear indicating a positive result. If the specimen does not contain Filariasis antibodies, a colored line will not appear in IgG/IgM region indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

Reagents

The test contains Filariasis antigen coated particles and anti-human IgG/IgM coated on the IgG/IgM region of membrane.

Precautions

- For laboratory *in vitro* diagnostic use. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against
- microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- The used test should be discarded according to local regulations.
- Humidity and temperature can adversely affect results.

Materials

Materials provided

- Test cassettes
- Dropers
- Buffer
- Product insert

Materials required but not provided

- Specimen collection containers
- Centrifuge
- Timer

Storage and Stability

Store as packaged in the sealed pouch at room temperature or refrigerated (35.6-86°F; 2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.


Collection and Storage of Specimens

- The Rapid Response™ Filariasis IgG/IgM Test Cassette can be performed using whole blood (from venipuncture), serum or plasma.
- Separate the serum or plasma from blood as soon as possible to avoid hemolysis. Only clear, non-hemolyzed specimens can be used.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 35.6-46.4°F (2-8°C) for up to 3 days.
- For long term storage, specimens should be kept below -4°F (-20°C). Whole blood collected by venipuncture should be stored at 35.6-46.4°F (2-8°C) if the test is to be run within 2 days of collection. Do not freeze whole blood specimens.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with federal regulations for transportation of etiologic agents.


Test Procedure

Allow the test, specimen, buffer and/or controls to reach room temperature (59-86°F; 15-30°C) prior to testing.


- Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
- Place the cassette on a clean and level surface. For **Serum or Plasma** specimen: Hold the dropper vertically and **transfer 1 drop of serum or plasma** (approximately 40µL) to the specimen area, then **add 2 drops of buffer** (approximately 80 µL), and start the timer, see illustration below. For **Venipuncture Whole Blood** specimen: Hold the dropper vertically and **transfer 1 drop of whole blood** (approximately 40µ) to the specimen area, then **add 2 drops of buffer** (approximately 80µL), and start the timer. See illustration below.
- Wait for the colored line(s) to appear. **Read the result at 15 minutes**, do not interpret the result after 20 minutes.




Results Interpretation



POSITIVE:* Two or three colored lines appear. One colored line should always appear in the control line region (C) and another one or two colored line(s) should be in the test line region(s) (IgM and/or IgG).




IgM Positive: A colored line appears in control region (C), another colored line appears in IgM region. It indicates a IgM positive test result for antibodies to Filariasis.




IgG Positive: A colored line appears in control region (C), another colored line appears in IgG region. It indicates a IgG positive test result for antibodies to Filariasis.


***NOTE:** The intensity of the color in the test line regions (IgM and IgG) may vary depending on the concentration of Filariasis antibodies present in the specimen. Therefore, any shade of color in the test line region (IgM and/or IgG) should be considered positive.



NEGATIVE: One colored line appears in the control region (C). No apparent colored line appears in the test region (T).



INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.



Quality Control

A procedural control is included in the test. A colored line appearing in the control region (C) is the internal procedural control. It confirms sufficient specimen volume and correct

procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

Limitations

- The Test Procedure and the Test Results Interpretation must be followed closely when testing the presence of antibodies to filarial parasites in serum, plasma or whole blood from individual subjects. Failure to follow the procedure may give inaccurate results.
- The Rapid Response™ Filariasis IgG/IgM Test Cassette is limited to the qualitative detection of antibodies to *W. bancrofti* and *B. malayi* in human serum, plasma or whole blood. The intensity of the test band does not have linear correlation with the antibody titer in the specimen.
- A negative result for an individual subject indicates absence of detectable *W. bancrofti* and *B. malayi* antibodies. However, a negative test result does not preclude the possibility of exposure to *W. bancrofti* and *B. malayi*.
- A negative result can occur if the quantity of *W. bancrofti* and *B. malayi* antibodies present in the specimen is below the detection limits of the assay, or the antibodies that are detected are not present during the stage of disease in which a sample is collected.
- Some specimens containing unusually high titer of heterophile antibodies or rheumatoid factor may affect expected results.
- The results obtained with this test should only be interpreted in conjunction with other diagnostic procedures and clinical findings.

Expected Values

The Rapid Response™ Filariasis IgG/IgM Test Cassette has been compared with Culture/Histology, demonstrating an overall accuracy of 97.0%.

Performance Characteristics

Clinical Sensitivity, Specificity and Accuracy

IgG Results

40 samples from patients with chronic lymphatic filariasis and 251 samples collected from a non-filariasis region were tested by the Rapid Response™ Filariasis IgG/IgM Test Cassette. Comparison for all subjects is shown in the following table:

Method	Results	ELISA		Total Results
		Positive	Negative	
Rapid Response™ Filariasis IgG/IgM Test Cassette	Positive	37	4	41
	Negative	3	247	250
Total Results		40	251	291

Relative Sensitivity: 92.5% (95%CI*: 79.6%-98.4%)

Relatively Specificity: 98.4% (95%CI*: 96.0%-99.6%)

Accuracy: 97.6% (95%CI*: 95.1%-99.0%)

*Confidence Interval

IgM Results

32 samples from patients with chronic lymphatic filariasis and 251 samples collected from a non-filariasis region were tested by the Rapid Response™ Filariasis IgG/IgM Test Cassette. Comparison for all subjects is shown in the following table:

Method	Results	ELISA		Total Results
		Positive	Negative	
Rapid Response™ Filariasis IgG/IgM Test Cassette	Positive	31	3	34
	Negative	1	248	249
Total Results		32	251	283

Relative Sensitivity: 96.9% (95%CI*: 83.8%-99.9%)

Relatively Specificity: 98.8% (95%CI*: 96.5%-99.8%)

Accuracy: 98.6% (95%CI*: 96.4%-99.6%)

*Confidence Interval

Precision

Intra-Assay

Within-run precision has been determined by using 10 replicates of four specimens: a negative, a low positive, a medium positive and a high positive. The negative, low positive, medium positive and high positive values were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the same four specimens: a negative, a low positive, a medium positive and a high positive. Three different lots of the Rapid Response™ Filariasis IgG/IgM Test Cassette have been tested using negative, low positive, medium positive and high positive specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity

Sera containing known amounts of antibodies to Filariasis have been tested with Hepatitis A, B, C, E, HIV and Syphilis. No cross-reactivity was observed, indicating that the Rapid Response™ Filariasis IgG/IgM Test Cassette has a high degree of specificity for antibodies to Filariasis.




Interfering Substances




The Rapid Response™ Filariasis IgG/IgM Test Cassette has been tested for possible interference from visibly hemolyzed and lipemic specimens, as well as serum specimens containing high bilirubin levels. In addition, no interference was observed in specimens containing up to 1,000 mg/dL hemoglobin; up to 1,000 mg/dL bilirubin; and up to 2,000 mg/dL human serum albumin.




Bibliography



- Lymphatic filariasis: the disease and its control. Fifth report of the WHO Expert Committee on Filariasis. WHO Tech Rep Ser 1992; 281-871.
- Michael E, Bundy DAP, Grenfell BT. Re-assessing the global prevalence and distribution of lymphatic filariasis. Parasitology 1996; 112:405-428.
- Eberhard ML, Lammie PJ. Laboratory diagnosis of filariasis. Clin. Lab Med 1991; 11:977-1010.
- More SJ, Copeman DB. A highly specific and sensitive monoclonal antibody-based ELISA for the detection of circulating antigen in bancroftian filariasis. Trop Med Parasitol 1990; 41:403-406
- Lammie PJ, Weil G, et al: Recombinant antigen-based antibody assays for the diagnosis surveillance of lymphatic filariasis-a multicenter trial. Flaria Jornal 2004: 3: 9-18.
- Baskar LK, Srikanth TR, et al: Development and evaluation of a rapid flow-through immunofiltration test using recombinant filarial antigen for diagnosis of brugian and bancroftian filariasis. Microbiol Immunol. 2004: 48: 519-25.

Glossary of Symbols

 Consult instructions for use
  Test per Kit
  Authorized Representative

 Store between 35.6°F to 86°F
  Use by
  Do Not Reuse

 Lot Number
  For *in vitro* diagnostic use only
  Catalogue #

 Do not use if package is damaged
  Manufacturer

 **BTNX Inc.**
 722 Rosebank Road,
 Pickering, ON L1W 4B2
 Canada



Technical Support: 1-888-339-9964