# Rapid Response®

## Fentanyl Home Test

REF FYL-1C156AHC

REF FYL-1C156AHC-1

### User Instructions


For Self-testing

A rapid test for the qualitative detection of FYL in human urine.

For self-testing *in vitro* diagnostic use.

### Materials Provided

### Materials required but not provided:

- Watch, clock or timer
- Urine collection cup

### Directions For Use

Allow the test to reach room temperature before testing. Wash your hands thoroughly.

#### 1 Collect Urine Sample

Collect urine in a clean and dry container.

#### 2 Testing

- Open the foil pouch, remove the test cassette and dropper, discard the desiccant.

- Place the test cassette on a clean and level surface.


- Hold the dropper vertically and transfer **3 full drops of urine** to the specimen well next to the letter "S" on the test cassette. **Do not spill urine into the test window.**

- Start a timer.

### Directions For Use

#### 3 Read Results

Read the test results at **5 minutes**, do not read results after 10 minutes.

### Interpretation of Results

#### NEGATIVE

- Two coloured lines appear.**
- One line appears next to the letter 'C' and the other line appears next to the letter 'T'.

#### POSITIVE

- Only one coloured line appears next to the letter 'C'. No coloured line appears next to the letter 'T'.**

#### INVALID

- Control line fails to appear.**

**IMPORTANT:** The result you obtained is called preliminary for a reason. The sample must be tested by laboratory in order to confirm the result.

### Intended Use

The Rapid Response® Fentanyl Home Test is a competitive binding, lateral flow immunochromatographic assay for qualitative detection of Fentanyl in human urine at the cutoff concentrations of 1 ng/mL. For self-testing *in vitro* diagnostic use. The test provides only preliminary test results. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. GC/MS or LC/MS is the preferred confirmatory method.

**WARNING: This testing device is not licensed for drug checking prior to consumption. When used for this purpose, there is a risk of inaccurate results.**

### Summary

Fentanyl, belongs to powerful narcotic analgesics, and is a  $\mu$  special opiates receptor stimulant. Fentanyl is one of the varieties that been listed in management of United Nations "Single Convention of narcotic drug in 1961". Among the opiates agents that under international control, fentanyl is one of the most commonly used to cure moderate to severe pain.<sup>1</sup> After continuous injection of fentanyl, the sufferer will have the performance of protracted opioid abstinence syndrome, such as ataxia and irritability etc,<sup>2,3</sup> which presents the addiction after taking fentanyl in a long time. Compared with drug addicts of amphetamine, drug addicts who take fentanyl mainly have got the possibility of higher infection rate of HIV, more dangerous injection behavior and more lifelong medication overdose.<sup>4</sup>

The Rapid Response® Fentanyl Home Test is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Fentanyl in urine. The Rapid Response® Fentanyl Home Test yields a positive result when Fentanyl in urine exceeds 1 ng/mL.

### Principle

The Rapid Response® Fentanyl Home Test is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Fentanyl, if present in the urine specimen below 1 ng/mL, will not saturate the binding sites of antibody-coated particles in the test device. The antibody-coated particles will then be captured by immobilized FYL conjugate and a visible coloured line will show up in the test line region. The coloured line will not form in the test line region if the FYL level exceeds 1 ng/mL because it will saturate all the binding sites of anti-FYL antibodies.

A drug-positive urine specimen will not generate a coloured line in the test line region, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a coloured line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

### Reagents

The test contains mouse monoclonal anti-FYL antibody-coupled particles and FYL-protein conjugate. A goat antibody is employed in the control line system.

### Precautions

- Read the entire user instruction prior to performing test.
- For self-testing *in vitro* diagnostic use only.
- For external use only.
- Do not use the test after expiration date printed on the package.
- Do not use the test if its foil pouch is torn or damaged.
- For single use. Discard after first use.
- The test device should remain in the sealed pouch until use.
- Contaminated or tainted urine sample may give false results.
- Keep out of reach of children.

### Storage and Stability

Store as packaged in the sealed pouch at 35.6-86°F (2-30°C) up to the date of expiration. Keep away from direct sunlight, moisture and heat. **DO NOT FREEZE.**

### Specimen Collection and Preparation

#### Urine Assay

- The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used.

#### Specimen Storage

- Urine specimens may be stored at 25°C/77°F for up to 24 hours or at 2-8°C/35.6-46.4°F for up to 48 hours prior to testing. For long-term storage, specimens may be frozen and stored at -20°C/-4°F for up to 12 months. Do not use the sample if it has been frozen and thawed more than 3 times.

### Performance Characteristics

#### Accuracy

About 80 clinical urine specimens fentanyl with known LC/MS values and tested by Rapid Response® Fentanyl Home Test. Each test was performed by three operators. Results were as follows:

Site	Device Result	Concentration by LC/MS (ng/mL)			
		≤ -50% Cut-off	-50% Cut-off to the cut-off	Cut-off to +50% cut-off	≥ +50% Cut-off
Site1	Positive	0	3	22	16
	Negative	25	12	2	0
	Total Results	100%	80.0%	91.7%	100%
Site2	Positive	0	3	22	16
	Negative	25	12	2	0
	Total Results	100%	80.0%	91.7%	100%

Site	Device Result	Concentration by LC/MS (ng/mL)			
		≤ -50% Cut-off	-50% Cut-off to the cut-off	Cut-off to +50% cut-off	≥ +50% Cut-off
Site3	Positive	0	2	21	16
	Negative	25	13	3	0
	Total Results	100%	86.7%	87.5%	100%

### Analytical Specificity

The following table list compounds that are positively detected in urine by Rapid Response® Fentanyl Home Test.

Fentanyl (Cut-off=1 ng/mL)	Concentration (ng/mL)
Acetyl fentanyl	1
Acrylfentanyl	1
ω-1-Hydroxyfentanyl	20,000
Isobutyryl fentanyl	1
Ocfentanil	2.5
Butyryl fentanyl	2.5
Furanyl fentanyl	5
Valeryl fentanyl	10
(±) β-hydroxythiofentanyl	2
4-Fluoro-isobutyrylfentanyl	50
Para-fluorobutyryl fentanyl	4
Para-fluoro fentanyl	3
(±)-3-cis-methyl fentanyl	50
Carfentanil	2
Sufentanil	7.5
Alfentanil	5,000
Despropionyl fentanyl (4-ANPP)	2,500
Trazodone	1,000
Remifentanil	>100 µg/mL
Norfentanyl	>100 µg/mL
Acetyl norfentanyl	>100 µg/mL
Norcarfentanil	>100 µg/mL

The following opioids compounds were tested at a concentration of 100 µg/mL. Negative results were obtained for all these compounds. There is no cross-reactivity for these compounds using the Rapid Response® Fentanyl Home Test.

6-Acetyl morphine	Naloxone
Amphetamine	Naltrexone
Buprenorphine	Norbuprenorphine
Buprenorphineglucuronide	Norcodeine
Codeine	Norketamine
Dextromethorphan	Normeperidine
Dihydrocodeine	Normorphine
EDDP	Noroxycodone
EMDP	Oxycodone
Fluoxetine	Oxymorphone
Heroin	Pentazocine (Talwin)
Hydrocodone	Pipamperone
Hydromorphone	Risperidone
Ketamine	Tapentadol
Levorphanol	Thioridazine
Meperidine	Tilidine
Methadone	Tramadol
Morphine	Tramadol-O- Desmethyl
Morphine-3-glucuronide	Tramadol-N- Desmethyl

## Precision

This study is performed by three POC personnel who don't know the sample number system participate in the study, at 3 POC site. Three lots were run in consecutive business days at each concentration for each lot per day. The results are as follows:

Concentration	n	Lot 1		Lot 2		Lot 3	
		+	-	+	-	+	-
0 ng/mL	60	0	60	0	60	0	60
0.25 ng/mL	60	0	60	0	60	0	60
0.5 ng/mL	60	0	60	0	60	0	60
0.75 ng/mL	60	2	58	1	59	0	60
1 ng/mL	60	35	25	28	32	26	34
1.25 ng/mL	60	60	0	60	0	60	0
1.5 ng/mL	60	60	0	60	0	60	0
1.75 ng/mL	60	60	0	60	0	60	0
2 ng/mL	60	60	0	60	0	60	0

## Effect of Urinary Specific Gravity

Total 12 urine samples of specific gravities (SG) ranging from 1.000-1.035 were collected. Values of SG levels were determined by a refractometer. Target drugs were spiked to these urine samples at +50% cut-off and -50% cut-off concentrations. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

## Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 4 to 9 in 1 pH unit increments and spiked with Fentanyl to 0.5 ng/mL and 1.5 ng/mL. The spiked, pH-adjusted urine was tested with the Rapid Response® Fentanyl Home Test in duplicate. The results demonstrated that varying ranges of pH do not interfere with the performance of the test.

## Interference

Potential interfering substances found in human urine of physiological or pathological conditions were added to drug-free urine and target drug fentanyl urine with concentrations at 50% below and 50% above cut-off levels. There urine samples were tested using three batches of each device. Compounds that showed no interference at a concentration of 100 µg/mL or specified concentrations are summarized in the following tables.

### Non Interference Compounds

Acetaminophen	Acetone (1000mg/dL)
Acetophenetidin	Acetylsalicylic acid
Albumin (100mg/dL)	Albuterol
Aminopyrine	Amitriptyline
Amobarbital	Amoxicillin
Ampicillin	Apomorphine
Ascorbic acid	Aspartame
Atropine	Benzilic acid
Benzoic acid	Benzoylcegonine
Bilirubin	Boric Acid (1%)
Bupropion	Caffeine
Carbamazepine	Chloral hydrate
Chloramphenicol	Chlorothiazide
Chlorpromazine	Cholesterol
Clomipramine	Clonidine
Cortisone	Cotinine
Potassium Chloride(1mg/mL)	Creatinine
Cyclobenzaprine	Deoxycorticosterone
Desipramine	Dextromethorphan
Diclofenac	Diflunisal
Digoxin	Diphenhydramine
DL-Tryptophan	DL-Tyrosine
Doxepin	Ecgonine methyl ester
Ephedrine	Erythromycin

Ethanol (1%)	Fenoprofen
Fluphenazine	Furosemide
Galactose (10mg/dL)	Gamma Globulin (500mg/dL)
Genitisc acid	Glucose (3000mg/dL)
Hemoglobin	Hydralazine
Hydrochlorothiazide	Hydrocortisone
Hydroxytyramine	Ibuprofen
Imipramine	Isoproterenol
Isosuxuprine (10µg/mL)	Sodium Chloride(1.6mg/mL)
Ketamine	Ketoprofen
Labeltalol	Lidocaine
Loperamide	Maprotiline
Meperidine (50µg/mL)	Meprobamate
Methapyrilene	Methaqualone
Methoxyphenamine	Metronidazole (300ug/ml)
N-Acetylprocainamide	NaCl (4000mg/dL)
Nalidixic acid	Naloxone
Naltrexone	Naproxen
Niacinamide	Nicotine
Nifedipine	Norethindrone
Nortriptyline	Noscapine
O-Hydroxyhippuric acid	Octopamine
Oxalic acid (100mg/dL)	Oxazepam
Oxolinic acid	Oxymetazoline
Papaverine	Penicillin G
Perphenazine	Phencyclidine
Phenelzine	Phenobarbital
Prednisone	Propoxyphene (50µg/mL)
Propranolol	Pseudoephedrine
Quinine	Ranitidine
Riboflavin (10mg/dL)	Salicylic acid
Secobarbital	Serotonin (5-Hydroxytyramine)
Sulfamethazine	Sulindac
Tetrahydrocortisone 3-(β-Dglucuronide)	Tetrahydrocortisone 3-acetate
Tetrahydrozoline	Thiamine
Thioridazine	Triamterene
Trifluoperazine	Trimethoprim
Tyramine	Urea (2000mg/dL)
Uric acid	Valproic acid (250µg/mL)
Venlafaxine	Verapamil
Zomepirac	β-Estradiol

## Q&A

### 1. How do I know if the Test worked well?

When the control line (C) appears, it means that the test unit is working well.

### 2. How soon can I read my results?

You can read your results after 5 minutes as long as a red line or pink coloured line has appeared next to the Control region (C), do not read results after 10 minutes.

### 3. When is the best time to run the test?

You can run the test at any time of the day.

### 4. How to read the test if the colour and the intensity of the lines are different?

The colour and intensity of the lines have no importance for result interpretation. The test should be considered as negative whatever the colour intensity of the test line (T) is.

### 5. What Is A False Positive Test?

A false positive test result means the drug is not present but shows detected by the device. The most common causes of a false positive test are cross reactants.

### 6. What Is A False Negative Test?

A false negative test means the drug is present but is not detected by the device. If the sample is diluted, or the sample is contaminated that may cause a false negative result.

## Limitations

1. The Test is for use with human urine only.
2. The Test only provides a preliminary result.
3. Substances, such as bleach and/or alum, in urine specimens may produce incorrect results. If adulteration is suspected, the test should be repeated with another urine specimen.
4. Test does not distinguish between drugs of abuse and certain medications.
5. A positive test result may be obtained from certain foods or food supplements.

## Bibliography

1. International Narcotics Control Board. Report of the International Narcotics Control Board for 2009[R]. New York: UN, 2010
2. Lane JC, Tension MB, Lawless ST, et al. Movement disorder after withdrawal of fentanyl infusion. J Pediatr, 1991, 119 (4): 649-651
3. Dominguez KD, Lomako DM, Katz RW, et al. Opioid withdraw in critically ill neonates. Ann Pharmacotherm, 2003, 37 (4): 473-477
4. European Monitoring Centre for Drugs and Drug Addiction. Annual Report 2009[R]. Lisbon: EMCDDA, 2010

## Glossary of Symbols

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Caution
	Temperature limit		Use by date		Manufacturer
	In vitro diagnostic medical device		Batch code		Do Not Reuse
	Do not use if package is damaged and consult instructions for use		Catalogue number		

**BTNX, Inc.**  
722 Rosebank Road  
Pickering, ON L1W 4B2 Canada  
Technical Support: 1-888-339-9964

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