 UTI-1SCH1

Rapid Response®

Urinary Tract Infection Screening Strip For Self Testing

Product Insert

For *in vitro* Diagnostic Use Only.

Unfold leaflet and flip this sheet over for instructions.

Intended Use

The Rapid Response® Urinary Tract Infection Screening Strip is a firm plastic strip onto which several separate reagent areas are affixed. The test is for the qualitative detection of the following analytes in urine: Leukocytes, Blood, Nitrite and Protein. The Rapid Response® Urinary Tract Infection Screening Strip is for single use in self-testing. The results should be interpreted in combination with other clinical information. Consult a physician before making any medical decision.

Summary

A urinary infection represents the most common disease of the urinary tract which includes the urethra, the bladder, the ureter and the kidneys. Men, women and children are likely to experience a urinary infection. It's mostly women who suffer from urinary infections, since the short urethra favours the penetration of germs. However, elderly males are also affected if they have an enlarged prostate which obstructs the urine flow. In healthy people, urine is sterile (i.e. it doesn't contain any micro-organisms). One of the best ways to keep your urinary tract sterile is to empty your bladder completely at regular intervals. Generally, an infection starts in the urethra and may then spread into the upper urinary tract as far as the kidneys. The symptoms vary considerably: burning when emptying the bladder, or a strong urge to urinate. The urine may also be cloudy or have a strong odour.

Principle

Leukocytes: This test reveals the presence of granulocyte esterases. The esterases cleave a derivatized pyrazole amino acid ester to liberate derivatized hydroxyl pyrazole. This pyrazole then reacts with a diazonium salt to produce a beige-pink to purple colour.

Blood: This test is based on the peroxidase-like activity of hemoglobin which catalyzes the reaction of diisopropylbenzene dihydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting colour ranges from orange to green to dark blue.

Nitrite: This test depends upon the conversion of nitrate to nitrite by the action of Gram negative bacteria in the urine. In an acidic medium, nitrite in the urine reacts with p-arsanilic acid to form a diazonium compound. The diazonium compound in turn couples with 1 N-(1-naphthyl) ethylenediamine to produce a pink colour.

Protein: This reaction is based on the phenomenon known as the "protein error" of pH indicators (Tetrabromophenol Blue), The anion produced by pH indicators under given conditions combines with the

cation produced by protein, then pH indicators Colours from yellow to green-blue for positive results.

Storage and Handling

Store as packaged at room temperature or refrigerated (2-30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

Warnings and Precautions

Read the instructions carefully before performing the test.

- For self-testing *in vitro* diagnostic use only.
- Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.
- A clean container - uncontaminated by cleaning fluids - to collect urine.
- Keep out of the reach of children.
- Do not use after the expiry date or if the pouch is damaged.
- Follow the indicated time strictly.
- Use the test only once. Do not dismantle and touch the reagent areas of the test strip.
- For external use only.
- The used test should be discarded according to local regulations.
- In case of difficulties in colour identification (such as Daltonism), ask for help in test reading.
- It is recommended to use the device and results interpretation under adequate lighting conditions; do not read the results in direct sunlight.

Technical Notes on Parameters

The test detects **LEUKOCYTES, BLOOD, NITRITE** and/or **PROTEINE** in urine.

- Leukocytes:** The presence of leukocytes in urine is an important symptom of an inflammation of the kidneys and the urinary tract, protein react with the pad and changing its colour to purple. When taking cephalixin and cephalothin, or high concentration of oxalic acid may also cause test results to be artificially low. Tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. High urinary protein may diminish the intensity of the reaction colour.
- Blood:** A uniform green colour conversion indicates the presence of hemoglobin or hemolyzite erythrocytes; scattered or compacted green spots indicate intact erythrocytes. General urine occult blood attributed to the following three reasons, one is a stone, inflammation, and cancer. On the aspect of inflammation, such as glomerulonephritis, pyelonephritis, cystitis, but may have hematuria, urine occult blood will have occurred. Whether stones kidney, ureter or bladder stones, may cause other situations such as occult blood. The tumor can also cause occult blood, such as benign or malignant tumor of kidney, ureter and bladder. Menstrual period, constipation may cause a positive result.
- Nitrite:** Gram-negative bacteria in urine convert nitrate from food into nitrite. Nitrite reacts with a chemical in the test field and leaves a pink shade. The test result may be distorted if urine does not stay for long in the bladder, due to hunger, a vegetable-free diet or antibiotic treatment. Comparing the test on a white background may aid in the detection of low nitrite levels, which might otherwise be missed.
- Protein:** An indicator on the test field reacts with protein in the urine, changing its colour to green. They may be found where there is inflammation of the bladder or prostate or bleeding in the urinary tract. Infusions containing polyvinylpyrrolidone may yield a false positive result. Chemical components in the test fields must be viewed as potentially dangerous substances, although they present no hazard provided that all test components are used in accordance with these instructions.

Control Procedure

To serve as good quality control, the instruction must be followed closely when performing the test. Failure to directions in insert may yield to inaccurate results.

Limitations

Note: The Rapid Response® Urinary Tract Infection Screening Strip may be affected by substances that cause abnormal urine colour such as drugs containing azo dyes (e.g. Pyridium®, AzoGantrisin®, AzoGantanol®), nitrofurantoin (Microdantin®, Furadantin®), and riboflavin.¹ The colour development on the test pad may be masked or a colour reaction may be produced that could be interpreted as false results.or elevated glucose in urine may cause test results to be low.

- Leukocytes:** The result should be read at 2 min to allow for complete colour development. The intensity of the colour that develops is proportional to the number of leukocytes present in the urine specimen. High specific gravity or elevated glucose concentrations (≥ 2,000 mg/dL) may cause test results to be artificially low. The presence of cephalixin, cephalothin, or high concentrations of oxalic acid may also cause test results to be artificially low. Tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. High urinary protein may diminish the intensity of the reaction colour. This test will not react with erythrocytes or bacteria common in urine.¹
- Blood:** A uniform green colour indicates the presence of myoglobin, hemoglobin or hemolyzed erythrocytes.¹ Scattered or compacted green spots indicate intact erythrocytes. To enhance accuracy, separate colour scales are provided for hemoglobin and for erythrocytes. Positive results with this test are often seen with urine from menstruating females. It has been reported that urine of high pH reduces sensitivity, while moderate to high concentration of ascorbic acid may inhibit colour formation. Microbial peroxidase, associated with urinary tract infection, may cause a false positive reaction.² The test is slightly more sensitive to free hemoglobin and myoglobin than to intact erythrocytes.
- Nitrite:** The test is specific for nitrite and will not react with any other substance normally excreted in urine. Any degree of uniform pink to red colour should be interpreted as a positive result, suggesting the presence of nitrite. Colour intensity is not proportional to the number of bacteria present in the urine specimen. Pink spots or pink edges should not be interpreted as a positive result. Comparing the reacted reagent area on a white background may aid in the detection of low nitrite levels, which might otherwise be missed. Ascorbic acid above 30 mg/dL may cause false negatives in urine containing less than 0.05 mg/dL nitrite ions. The sensitivity of this test is reduced for urine specimens with highly buffered alkaline urine or with high specific gravity. A negative result does not at any time preclude the possibility of bacteruria. Negative results may occur in urinary tract infections from organisms that do not contain reductase to convert nitrate to nitrite; when urine has not been retained in the bladder for a sufficient length of time (at least 4 hours) for reduction of nitrate to nitrite to occur; when receiving antibiotic therapy or when dietary nitrate is absent.
- Protein:** Any green colour indicates the presence of protein in the urine. This test is highly sensitive for albumin, and less sensitive to hemoglobin, globulin and mucoprotein.¹ A negative result does not rule out the presence of these other proteins.

False positive results may be obtained with highly buffered or alkaline urine. Contamination of urine specimens with quaternary ammonium compounds or skin cleansers containing chlorhexidine may produce false positive results.¹ The urine specimens with high specific gravity may give false negative results.

Performance Characteristics

Parameters of importance to the user are sensitivity, specificity, accuracy and precision. Generally, this test has been developed to be specific for

the parameters to be measured with the exceptions of the interferences listed. Please refer to the Limitations section in this package insert. Interpretation of visual results is dependent on several factors: the variability of colour perception, the presence or absence of inhibitory factors, and the lighting conditions when the strip is read. Each colour block on the chart corresponds to a range of analyte concentrations.

Analytical Sensitivity:

Blood	5 Ery/μL
Protein	12 mg/dL
Nitrite	0.05 mg/dL
Leukocytes	10 Leu/μL

Precision: There is no significant difference within runs, between days, sites, lots and operators.

Accuracy: The difference between the test result and the labeled value of the corresponding reference solution is no more than 1 order of magnitude in the same direction and there is no reverse difference. Positive reference solution does not produce negative result, and negative reference solution does not produce positive result.

FAQ

What should I do if my test result is positive?

- Remember that a positive result doesn't mean that all four substances have been detected in your urine. Even if your result is positive for just one of them, it is most likely that something is wrong in your urine, even if the reason may not be a urinary infection. Get in touch promptly with your own doctor, who will be able to give a more accurate diagnosis. When you visit your doctor, please take these instructions with you so that he/she will be better informed as to the type of test you have performed.

What should I do if my test result is negative?

- Remember that your test result is only negative if the result on the test field for all four substances is negative. But if you still feel the signs of a UTI or have any other symptom, then contact your own doctor to arrange a more thorough examination.

Bibliography

- Henry JB, et al. Clinical Diagnosis and Management by Laboratory Methods, 20th Ed. Philadelphia. Saunders. 371-372, 375, 379, 382, 385, 2001.
- Ma Junlong, Cong Yulong. The effect of bacteriuria on the determination of urine red blood cells by urine analyzer. Chinese Journal of Medical Examination, 1999, 22(4): 205.
- Shuai Lihua, Jiujiang Medical Journal 2002, 17 (2): 122.

Glossary of Symbols

 Consult instructions for use

 Test per Kit

 Catalogue number

 Store between 2°C to 30°C

 Use by date

 Do Not Reuse

 *In vitro* diagnostic medical device

 Lot Number

 Manufacturer

 Do not use if package is damaged and consult instructions for use

BTNX Inc.
722 Rosebank Road,
Pickering, ON L1W 4B2 Canada
Technical Support: 1-888-339-9964
www.btnx.com
Monday-Friday 9:00AM - 5:30PM EST

Rapid Response®

Urinary Tract Infection Screening Strip

For Self Testing. For *in vitro* diagnostic use only.

Allow the strip and urine specimen to reach room temperature 15-30°C (59-86°F) prior to testing.

Notes

It's recommended to take a sample of urine for the test in the early morning since it's the most concentrated. The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances.

For women only: The test should not be performed during or for three days after your menstrual period. The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

Do not make any important medical decisions without first referring to your doctor.

Materials

Unpack the test components from the kit and make sure that all the packaging is intact.


For each test you will need:

Sealed Pouch
Containing Test Strip

Product Insert
(This Document)

Colour Chart

Plastic Cup


Timer or Watch with
Second Hand
(Not provided)

Specimen Container
(Not provided)

Instructions

Collect urine: It is recommended to collect midstream urine sample in the supplied plastic cup or using a clean cup without any residual detergents. Make sure to fill up the cup with urine.


Open the foil pouch and take out the test strip. **Do not touch the test fields.** Once opened the pouch, it is recommended to perform the test **immediately**.


Dip the test strip in the urine sample.


Note: Press the strip and make sure that all four test fields are immersed for about 1-2 seconds.

Then remove the test strip and wipe off any surplus urine against the rim of the container or with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas. Failure in wiping off any surplus urine may give false results. The strip should be maintained in the horizontal position facing upwards.

Wait for 2 min (do not read results after three minutes). Read the result separately for each parameter, compare colour with colour chart provided.


Reading the Results

Read the result separately for each parameter; compare colour with colour chart provided. Colour changes on the edges of the test pads or colour changes after more than 3 minutes have to be ignored.


NEGATIVE

The Test field for **leukocytes** stayed whitish.
The Test field for **blood** stayed mustard yellow.
The Test field for **nitrite** stayed white.
The Test field for **protein** stayed yellowish.


POSITIVE RESULT FOR LEUKOCYTES

If the colour of the test field has changed to **purple**, then leukocytes have been found in your urine.


POSITIVE RESULT FOR BLOOD

If the colour of the test field has changed to **green** (or some green spots appear on the background), then blood has been found in your urine.

POSITIVE RESULT FOR NITRITE

If the colour of the test field has changed to **pink**, then nitrites have been found in your urine.

POSITIVE RESULT FOR PROTEIN

If the colour of the test field has changed to **green**, then proteins have been found in your urine.

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