# Rapid Response™

## Urinary Tract Infection Screen Strips

REF UTI-1SHU1

REF UTI-1SHU3

### Product Insert


It is intended for over-the-counter home use only.  
For *in vitro* diagnostic use only.  
For urine testing only.

### Materials

#### Materials provided:

Test Strip

Product Insert

**Materials required but not provided:** Timer or a watch with a second hand.


### Procedure

**Before testing: It is recommended to take a sample of urine for the test in the early morning since it's the most concentrated.** The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances. **Do not make any important medical decisions without first referring to your doctor.**

**1** Remove the test strip from the foil pouch. Retain the foil pouch to use the color chart on the outside to interpret the results. Do not use the test strip if the pouch is torn or damaged.


**2** Hold the end of the test strip and begin urinating. After 1 or 2 seconds of passing urine, hold the strip downward with the test pads pointing into the urine stream for **1-2 seconds**. Make sure the test pads are completely wet.


### Procedure

**3** Immediately bring the edge of the test strip into contact with an absorbent material (e.g. paper towel) to remove excess urine. Lay the test strip with the test pads facing upwards and begin timing.

**4** Interpret results at 1 minute, do not read results after 3 minutes. Compare the test pads to the color blocks on the foil pouch. Hold the test strip close to the color blocks and carefully match each test pad to the color chart for that test.


### Interpreting Your Results

Interpret the results by comparing the test pads to the color blocks on the foil pouch. Match the color of the test pad to the closest color block on the foil pouch.

Test Strip Results		Recommendations
Leukocytes (LEU)	Nitrite (NIT)	
Negative (Neg.)	Negative (Neg.)	No sign of UTI detected. However, if you still have symptoms, consult your physician; there are some cases of UTI that the test may not detect.
Trace	Negative (Neg.)	Results suggest additional testing is required. Test again the following day with a new strip. If you have symptoms or if you get another trace leukocyte result, consult with your physician.
Negative or Trace	Positive	Results suggest a sign of a UTI. Consult your physician immediately.
Positive +, ++, +++	Positive	Results suggest a sign of a UTI. Consult your physician immediately.
Positive +, ++, +++	Negative (Neg.)	Results suggest a sign of a UTI. Repeat the test next time you urinate using a new test strip. Make sure to wash the genital area first. If LEU is still positive, consult your physician.

### Intended Use

The Rapid Response® Urinary Tract Infection Screen Strips is for the qualitative detection of Leukocytes and Nitrite in human urine as an aid in the screening of urinary tract infection (UTI). It is intended for over-the-counter home use only.

### Reagents

Leukocytes (LEU)	0.5% w/w pyrrole amino acid ester 0.4% w/w diazonium salt 99.1% w/w buffer and non-reactive ingredients
------------------	---

### Nitrite (NIT)

4.5% w/w p-aminobenzene sulfonamide  
1.5% w/w 1,2,3,4-tetrahydro-benzo(h)quinolinolone-3-pheno;  
94% w/w buffer and non-reactive ingredients

### Detection Time/Range

Item	Sensitivity	Reportable Range	Read Time
Nitrite (mg/dL)	0.05	Negative-Positive	1 minute
Leukocytes (Cells/μL)	15	Negative-500	1 minute

### Warnings and Precautions

**Read the instructions carefully before performing the test.**

- For *in vitro* diagnostic use only.
- For urine testing only. Do not use for blood testing.
- Keep out of the reach of children.
- Do not use after the expiry date or if the pouch is damaged.
- Follow the indicated time strictly.
- Use the test immediately or within 2 hours after opening.
- Use the test only once. Do not take apart or touch the test pad areas of the test strip.
- For external use only. Not to be taken internally.
- The used test should be discarded according to local regulations.
- Individuals with color-impaired vision may not be able to adequately interpret test results.
- Read the results in bright, natural light. Do not read under colored light sources to avoid color deviation.
- Any unused test strip with discoloration or darkening of the test pads should be considered damaged or defective and should be discarded.

### Storage and Stability

- Store test strips as packaged in the original sealed pouch.
- Store in a dry place at temperatures between 35.6-86°F (2-30°C) up to the date of expiration.
- Do not freeze.
- Keep away from direct sunlight.
- Strips can be used only once.

### Limitations

There is the possibility that this test may produce false results. The test may be affected by medication or dietary supplement that may change urine color. Do not use test strips if you observe abnormal urine colors (such as bright yellow or green). Do not use test strips while taking antibiotics. Consult with your physician before making any medical decision.

#### NOTE:

- Leukocytes: High glucose ( $\geq 1,000$  mg/dL) or high Vitamin C ( $\geq 150$  mg/dL) in urine may cause test results to be low. Urine with a  $\text{pH} \geq 8.0$  may cause a false positive leukocyte reading. Dehydration may cause false negative leukocyte reading.
- Nitrite: High Vitamin C ( $\geq 150$  mg/dL) in urine may cause a false negative nitrite reading. Any degree of uniform pink color should be read as a positive result. You may get a negative result if you have UTI caused by bacteria that do not change nitrate to nitrite, when urine has not been held in the bladder for more than 4 hours, when taking antibiotics, or when your diet does not include nitrates.

### Frequently Asked Questions

#### What are the common symptoms of UTI?

The common symptoms of UTI include burning or pain during urination, frequent urination, fever, and back or groin pain.

#### Can I have UTI without symptoms?

Yes. UTI can occur without symptoms in young children, pregnant women and the elderly.


#### Do I have to test with first morning urine?

Although you can test at any time of the day, your first morning urine would have the most bacteria in it if UTI is present. Avoid testing during your menstrual period.

#### Do certain foods, alcohol or common medications affect the test?

In general, most substances will not interfere with the test. Consult your physician if you are taking antibiotics. Some medications (e.g. Vitamin B2) can affect the color of your urine and may cause false results. Large quantities of Vitamin C may cause false negative results.

### Index of Symbols

**Manufactured for:**  
**BTNX Inc.**  
722 Rosebank Road,  
Pickering, ON L1W 4B2 Canada  
**Technical Support: 1-888-339-9964**  
**Email: support@btnx.com**  
**Mon.-Fri. 9:00AM-5:30PM EST**  
**www.btnx.com**

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