



Rapid Response® Plastic Collection Cups



Key Features

- Leak proof screw-on cap
- Made with polypropylene for clarity and durability
- Available with and without a temperature strip (ideal for drug screening)
- 120mL Capacity

Ideal For

- Use with Rapid Response® Drug Screening Strips, Panels, and Cassettes
- Use with Rapid Response® Urinalysis Reagent Strips
- Use with Rapid Response® HCG Pregnancy Test Strips and Cassettes

Ordering Information

| Cat. No. | Product Name | Description | Qty /Kit |
|-----------|---|--|----------|
| 02SPC-14 | Plastic Collection Cup | Plastic Sample Collection Cup with Cap | 25 |
| 02SPC-14A | Plastic Collection Cup With Temperature Strip | Plastic Sample Collection Cup with Cap and Temperature Strip | 25 |