

REF BZO-18S3-100


Rapid Response™ Benzodiazepine Test Strip

(Liquid / Powder)

Benzodiazepine are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Risk of physical dependence increases if Benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception. The detection period for the Benzodiazepines in the urine is 3-7 days.

The Rapid Response™ Benzodiazepine Test Strip (Liquid / Powder) is an easy-to-use, rapid lateral flow assay for the qualitative detection of benzodiazepines in liquid or powder samples.


Testing Procedure





For illustrative purposes only. Refer to product insert for more details.

Results Interpretation


Positive
Contains
BZO



Negative
BZO could not
be detected



Invalid



Ordering Information

Product Code	Product Name	Contents
BZO-18S3-100	Rapid Response™ Benzodiazepine Test Strip (Liquid / Powder)	100 Tests

To learn more contact your local **Sales Representative**, call us at **+1 888 506-2658**, or email us at **info@lochnessmedical.com**

Kit Contents

- Individually packaged test strips
- Product insert
- Results interpretation card

Product Details

- Detects Benzodiazepines at levels above 300 ng/mL (typical drug samples contain over 1000 times this cut-off)
- Results interpretation guide printed on each pouch
- Built-in procedural control
- Results in minutes
- Compact format for easy distribution

Detects Benzodiazepine Analogues including:

Bromazepam	Diazepam	Estazolam
------------	----------	-----------




To learn more about Rapid Response™ Benzodiazepine Test Strip (Liquid / Powder), scan this QR code.

Our Harm Reduction Line Up



Rapid Response™ Xylazine Test Strips

Product Code	Product Name	Contents
XYL-18S2-100	Rapid Response™ Xylazine Test Strip (Liquid / Powder)	100 Tests
XYL-18S2-10	Rapid Response™ Xylazine Test Strip (Liquid / Powder)	10 Tests
XYL-18S2-5K	Rapid Response™ Xylazine Test Strip Kit (Liquid / Powder)	5 Tests



Rapid Response™ Fentanyl Test Strips

Product Code	Product Name	Contents
FYL-18S7-100	Rapid Response™ Fentanyl Test Strip (Liquid / Powder)	100 Tests
FYL-18S7-10	Rapid Response™ Fentanyl Test Strip (Liquid / Powder)	10 Tests
FYL-18S7-5	Rapid Response™ Fentanyl Test Strip Kit (Liquid / Powder)	5 Tests



Rapid Response™ Benzodiazepine Test Strips

Product Code	Product Name	Contents
BZO-18S3-100	Rapid Response™ Benzodiazepine Test Strip (Liquid / Powder)	100 Tests



Rapid Response™ Nitazene Test Strips

Product Code	Product Name	Contents
NTZ-18S26-100	Rapid Response™ Nitazene Test Strip (Liquid / Powder)	100 Tests



Visit Our Harm Reduction Page

www.lochnessmedical.com/category/Harm-Reduction

To learn more contact your local **Sales Representative**, call us at **+1 888 506-2658**, or email us at info@lochnessmedical.com