

**DRUG SCREENING**

# Rapid Response™ Drug Screen Flat Cup

with a **Fentanyl Test Strip**

| For *in vitro* diagnostic use only.

- ✓ Multi-drug urine detection
- ✓ Fentanyl included (1 ng/mL cut-off)
- ✓ Results in 5 minutes
- ✓ Compatible with digital workflow solutions
- ✓ Built-in temperature strip and adulteration strips *(included in SKU: D12.49-1VAF)*



**FDA 510(k)  
Cleared &  
CLIA Waived**

for professional /  
healthcare use.



SKU	PRODUCT NAME	DRUG PARAMETERS
D12.49-1VAF	<b>12 Drug Screen Flat Cup</b> with Fentanyl + Adulteration Strips <i>CLIA WAIVED</i>	OXY100 BZO300 AMP500 COC150 THC50 MDMA500 BUP10 BAR300 MTD300 MET500 MOR300 <b>FYL1 + OXI/CRE/PH</b>
D13.5-1VF	<b>13 Drug Screen Flat Cup</b> with Fentanyl <i>CLIA WAIVED</i>	OXY100 BZO300 AMP500 COC150 THC50 MDMA500 BUP10 BAR300 MTD300 MET500 PCP25 OPI300 <b>FYL1</b>

## Workflow Solutions

### The RapidReader App



Scan QR code to learn more.



The RapidReader App helps treatment centers streamline their drug screening process through a compact and user-friendly mobile application. Using a cellphone or tablet, the app utilizes AI to capture individual test results from a photo and automatically uploads the report to the patient's profile in the EMR.

### Lateral Flow Scanner



Scan QR code to learn more.



The Rapid Response™ Lateral Flow Scanner is an ideal solution for treatment centers looking to streamline medium- to high-volume drug screen reporting. The system can simultaneously scan up to three patients' Rapid Response™ Drug Tests. Powered by AI, it captures the results and uploads a report directly to each patient's profile in the EMR system.



[bit.ly/4pW5gnr](https://bit.ly/4pW5gnr)

**SCAN QR CODE** to visit our website for more resources.  
 Contact us for pricing or bulk discounts at [orders@lochnessmedical.com](mailto:orders@lochnessmedical.com)  
 or call us at +1 888-506-2658



## Rapid Response™ Drug Screen Flat Cup

The Rapid Response™ Drug Screen Flat Cup is a CLIA Waived, one-step rapid test that detects multiple drugs and metabolites in urine and is now enhanced with the addition of a fentanyl test strip. With the rising prevalence of fentanyl in the drug supply, incorporating fentanyl testing into our CLIA Waived test cups is a critical feature, providing you with accurate and rapid results that meet your testing needs.

**Sample:** Urine

**Format:** Cup

**Time-to-Result:** 5 minutes

**Storage Condition:** 35.6-86°F / 2-30°C

**Test Principle:** Lateral Flow Assay

### Detects Multiple Drugs

This CLIA-Waived drug screening cup detects a wide range of drugs in various combinations, including fentanyl. It's an ideal choice for professionals seeking a precise and reliable drug screening solution.

### Point-of-Care Testing

Screen for multiple drugs and metabolites in patient urine samples in just 5 minutes.

### Easy & Convenient

A simple one-step test with integrated drug test strips enables easy collection and testing. Each cup features a leak-proof cap and a temperature strip for added reliability.

### Kit Content

Each kit includes test cups, an adulteration color chart (when applicable), and a product insert.

## WHY CHOOSE RAPID RESPONSE™



### Fentanyl Testing Integrated

Fentanyl is a potent opioid increasingly found in the drug supply, making rapid detection critical. These test cups include a highly sensitive fentanyl test strip, with a low cut-off of 1 ng/mL, allowing quick testing without the need to test for fentanyl separately.



### CLIA Waived

These test cups are FDA 510(k) cleared for professional or other healthcare use under a CLIA Certificate of Waiver for simplified testing.



### Built-in Tamper / Validity Checks

Verify specimen quality with the built-in temperature strip and adulteration strips (included in SKU: D12.49-1VAF) to prevent against sample adulteration and tampered results.



### Rapid, One-Step Workflow

Collect and test multiple drugs in a single cup, no extra strips or accessories needed, enabling fast, on-site clinical decisions with results in just 5 minutes.



### Digital Tools Compatible

The flat cup format design pairs well with the Lateral Flow Scanner and The RapidReader App to help boost efficiency and improve record-keeping.



### High Quality Standards

Rapid Response™ is a trusted brand that is backed by rigorous regulatory validation, independent studies, ongoing proficiency testing, and maintains the highest internationally quality standards audited annually by Intertek.