



Distributed by



DASH® Rapid PCR System for Urgent Care

Rapid PCR Made Simple™



Lab-Quality Results in 15 Minutes at the Point of Care

- ✓ **SIMPLE**
Intuitive, easy-to-operate user experience
- ✓ **RAPID**
Delivers results in 15 minutes, with only one minute of hands-on time
- ✓ **AFFORDABLE**
Low-cost solution enables PCR testing at the point of care
- ✓ **ACCURATE**
Lab-quality results enhance confidence in treatment decisions

Test and Treat in a Single Patient Visit



PATIENT ARRIVES AND
SAMPLE IS COLLECTED



RUN DASH® SARS-CoV-2
& FLU A/B TEST



RESULTS IN JUST
15 MINUTES



CLINICIAN PRESENTS RESULTS
AND DISCUSSES TREATMENT PLAN

Full Reimbursement Potential with Existing CPT Code

The DASH® SARS-CoV-2 & Flu A/B Test is affordably priced to support high-quality PCR testing at the point of care. The test utilizes the existing CPT code 87636-QW.





Nuclein

Distributed by



LOCHNESS
MEDICAL

DASH® SARS-CoV-2 & Flu A/B Test

FDA Cleared & CLIA Waived | 15-Minute Runtime | PCR at the POC

Results From a Prospective Clinical Study Against an FDA-Cleared PCR Test

POSITIVE AGREEMENT

(Sensitivity)

SARS-CoV-2	95.2%
Flu A	94.3%
Flu B	97.3%

NEGATIVE AGREEMENT


(Specificity)

SARS-CoV-2	99.5%
Flu A	98.1%
Flu B	99.2%

Source: DASH® SARS-CoV-2 & Flu A/B Test package insert.

Why Choose the DASH® Rapid PCR System?

- ✓ Rapid PCR test results at the point of care
- ✓ Lab-quality sensitivity and specificity
- ✓ Anterior nasal swab specimens
- ✓ Qualitative, intuitive patient report
- ✓ Wi-Fi and cloud connectivity enabled



DASH® Rapid PCR System Capabilities



Cartridges are stored at room temperature



Designed to support a wide variety of sample types



Intuitive operation by non-laboratory personnel



No calibration or preventive maintenance required



Expansive planned test menu, including both viral and bacterial pathogens



Built to be multiplex-capable in a single cartridge

orders@lochnessmedical.com | 1-888-506-2658 | www.lochnessmedical.com