

REF DX.X-2FC1

Rapid Response™

Multi-Drug Saliva Test Cup

(Saliva)

DRUG SCREENING

The Rapid Response® Saliva Test Cup is a rapid, non-invasive drug screening solution that can simultaneously detect multiple drugs in saliva samples. The test is available in a variety of drug parameter and cut-off combinations and can test up to 16 parameters in a single cup. The Saliva Test Cup is also compatible with our RapidReader App allowing you to upload results to your ATS/EMR/EHR system in seconds.



Product Details

The Saliva cup is a rapid non-invasive all-in-one collection swab and drug test solution that can simultaneously detect multiple drugs in saliva.

- No contact with urine
- Easy to monitor collection to prevent adulteration
- Difficult to falsify test results
- Results in 10 minutes

Sample: Saliva


Format: Cup

Time-to-Result: 10 minutes

Storage Condition: 36-86°F/2-30°C

Test Principle: Lateral Flow Assay

Detection Window: Up to 24 hours before collection



Now compatible with our

RapidReader App

Streamline your drug screening process by alleviating the burden of physical paperwork and manual inputs. The RapidReader App is a reliable solution to support users with results collection, reporting, and upload directly from any mobile device.



Scan to see the app in action


Multi-Drug Saliva Test Cup



Easy to Use

- 1 Use the collection swab to sweep inside of mouth until colour appears on the saturation indicator strip.
- 2 Place the collection swab into the test cup and screw the cap back on.
- 3 Read results at 10 minutes.

For illustrative purposes only, please read the product insert carefully before use.



Popular Parameter Combinations

Each Rapid Response® Saliva Test Cup product has a set combination of parameters represented by a drug parameter code. The first three letters of a drug parameter code refer to the parameter detected; the number that follows represents the Cut-off Level in ng/mL for that parameter.

REF No.	Product Details	Drug Parameters	Quantity
D5.1.1-2FC1	5 Parameter Saliva Test Cup	THC40, COC20, AMP50, MET50, OPI40	25 per Kit
D6.3.1-2FC1	6 Parameter Saliva Test Cup	COC20, MET50, PCP10, THC40, OPI40, AMP50	25 per Kit
D6.3-2FC1	6 Parameter Saliva Test Cup	COC20, MET50, PCP10, THC30, OPI40, AMP50	25 per Kit
D10.3.1-2FC1	10 Parameter Saliva Test Cup	THC40, COC20, OPI40, AMP50, MET50, BAR50, BZO10, MDMA50, MTD30, PCP10	25 per Kit
D10.32.1-2FC1	10 Parameter Saliva Test Cup	COC20, OPI40, AMP50, MET50, BAR50, BZO10, MDMA50, MTD30, PCP10, OXY20	25 per Kit
D12.5.1-2FC1	12 Parameter Saliva Test Cup	THC40, COC20, OPI40, AMP50, MET50, BAR50, BZO10, MDMA50, MTD30, PCP10, BUP5, OXY20	25 per Kit

Customizable parameter and cut-off combinations are available for qualifying orders.

To learn more contact your local Sales Representative, call us at +1-888-339-9964, or email us at sales@btnx.com.