

**REF XYL-18S2-5K**

# Rapid Response™ Xylazine Test Strip Kit


(Liquid / Powder)

New enhanced kit. Same reliable strip.  
—Testing just got easier!

Xylazine is a non-opioid tranquilizer used as a sedative, analgesic, and muscle relaxant in veterinary medicine. Recently, it has commonly been encountered in the recreational drug supply, most commonly in opioid “dope” (e.g. fentanyl, heroin). Combining xylazine with other drugs that cause central nervous system depression compounds the sedative effects and can increase the risk of overdose and death.


The Rapid Response™ Xylazine Test Strip Kit (Liquid / Powder) is an easy to use, rapid lateral flow assay for the qualitative detection of xylazine in liquid or powder dope samples. For forensic use only. Not an IVD.

## Testing Procedure



For illustrative purposes only. Refer to product insert for more details.

## Results Interpretation



## Kit Contents

- 5 Individually packaged test strips
- 5 Buffer tubes
- 5 Swabs
- Tube stand
- Results interpretation card
- Instructions for use

## Product Details

- Detects Xylazine at concentrations above 1000 ng/mL
- Results interpretation guide printed on the back of each pouch
- Includes all materials required to test for Xylazine in drug samples



bit.ly/3G5RmuX

## Read the CFSRE Report

Rapid Response™ Xylazine Test Strip (Liquid / Powder)'s performance was deemed acceptable for drug checking purposes, demonstrating high sensitivity, specificity, and precision when evaluated by the Center for Forensic Science Research and Education (CFSRE) and the Philadelphia Department of Public Health (PDPH).



www.btnx.com/xylazine-kit

## Visit Our Website for more Harm Reduction Resources

To learn more about Rapid Response™ Xylazine Test Strip (Liquid / Powder), scan this QR code.

## Ordering Information

Product Code	Product Name	Contents
XYL-18S2-5K	Rapid Response™ Xylazine Test Strip Kit (Liquid / Powder)	5 Tests


To learn more contact your local **Sales Representative**, call us at **+1 888 506-2658**, or email us at **info@lochnessmedical.com**

## Our Harm Reduction Line Up



### Rapid Response™ Xylazine Test Strips

Product Code	Product Name	Contents
XYL-18S2-100	Rapid Response™ Xylazine Test Strip (Liquid / Powder)	100 Tests
XYL-18S2-10	Rapid Response™ Xylazine Test Strip (Liquid / Powder)	10 Tests
XYL-18S2-5K	Rapid Response™ Xylazine Test Strip Kit (Liquid / Powder)	5 Tests



### Rapid Response™ Fentanyl Test Strips

Product Code	Product Name	Contents
FYL-18S7-100	Rapid Response™ Fentanyl Test Strip (Liquid / Powder)	100 Tests
FYL-18S7-10	Rapid Response™ Fentanyl Test Strip (Liquid / Powder)	10 Tests
FYL-18S7-5	Rapid Response™ Fentanyl Test Strip Kit (Liquid / Powder)	5 Tests



### Rapid Response™ Benzodiazepine Test Strips

Product Code	Product Name	Contents
BZO-18S3-100	Rapid Response™ Benzodiazepine Test Strip (Liquid / Powder)	100 Tests



### Rapid Response™ Nitazene Test Strips

Product Code	Product Name	Contents
NTZ-18S26-100	Rapid Response™ Nitazene Test Strip (Liquid / Powder)	100 Tests



**Visit the BTNX Harm  
Reduction Resource Page**  
[www.btnx.com/HarmReduction](http://www.btnx.com/HarmReduction)

To learn more contact your local **Sales Representative**, call us at **+1 888 506-2658**, or email us at [info@lochnessmedical.com](mailto:info@lochnessmedical.com)